BÜHLMANN GanglioCombi™ *Light* ELISA

with enzyme labels IgG/IgM Mix, IgG and IgM

Detection of anti-ganglioside
auto-antibodies by ELISA
(GM1, GD1b, and GQ1b)

**For research use only.
Not for use in diagnostic procedures.**

EK-GCL-S-U 96 wells

Revision date: 2017-06-09

ENGLISH

INTENDED USE

BÜHLMANN GanglioCombi™ *Light* ELISA, is an *in vitro* test intended to detect antibodies against defined relevant neural antigens GM1, GD1b and GQ1b (IgG and IgM) in serum samples.

This product is for research use only. It is not intended for use in diagnostic procedures.

INTENDED APPLICATION

With regard to the 3 different enzyme labels, the device components allow three application options:

1. Testing with the IgG/IgM mix conjugate allows to screen for the presence of antibodies.
2. Testing with individual IgG and/or IgM conjugates for antibody isotype determination.
3. For laboratory work-up we suggest combining both procedures: sample screening using the mix (option 1), followed by differentiation of mix-positive samples by separate IgG and IgM enzyme labels (option 2).

PRINCIPLE OF THE ASSAY

BÜHLMANN GanglioCombi™ *Light* ELISA is based on the enzyme-immunometric assay technique. The wells of the provided microtiter plate are coated with gangliosides: GM1, GD1b and GQ1b.

Calibrator, controls, and sera are incubated in the microtiter wells and anti-ganglioside antibodies present in the samples bind to the immobilized gangliosides. After washing off unbound substances, the antibodies are detected with horseradish-peroxidase (HRP) labelled antibodies against human IgG and/or IgM. Following a second washing step in which unbound enzyme label is removed, a substrate solution containing tetramethyl-benzidine (TMB) is added. A blue colour develops in proportion to the amount of antibodies bound to the immobilized gangliosides. Colour development is stopped by adding an acidic stop solution (diluted sulphuric acid) which turns the blue solution into yellow. The intensity of the colour is measured at 450 nm.

The measured absorbance is proportional to the titre of antibodies present in a given sample. The titres of antibodies are expressed as % ratios of the calibrator.

REAGENTS SUPPLIED AND PREPARATION

Reagents	Quantity	Code	Reconstitution
Microtiter Plate precoated with gangliosides	12 x 8 wells	B-GCL-MP	Ready to use
Plate Sealer	3 pieces		
Wash Buffer Concentrate (10X) with preservatives	1 bottle 100 mL	B-GCO-WB	Dilute with 900 mL of deionized water
Incubation Buffer with preservatives	1 bottle 100 mL	B-GCO-IB	Ready to use
Calibrator Lyophilized with preservatives	1 vial	B-GCO-CA	Add 1.5 mL of Incubation Buffer
Negative, Low and Medium Control Lyophilized with preservatives	3 vials	B-GCO-CONSET	Add 1.5 mL of Incubation Buffer

Reagents	Quantity	Code	Reconstitution
Enzyme Label IgG Anti-human IgG Ab conjugated to HRP in a protein-based buffer with preservatives	1 vial 11 mL	B-GCO-ELG	Ready to use
Enzyme Label IgM Anti-human IgM Ab conjugated to HRP in a protein-based buffer with preservatives	1 vial 11 mL	B-GCO-ELM	Ready to use
Enzyme Label Mix Anti-human IgG and IgM Ab conjugated to HRP in a protein-based buffer with preservatives	1 vial 11 mL	B-GCO-ELGM	Ready to use
TMB Substrate TMB in citrate buffer	1 vial 11 mL	B-TMB	Ready to use
Stop Solution 0.25 M sulfuric acid	1 vial 11 mL	B-STTS	Ready to use Corrosive agent

Table 1

STORAGE AND SHELF LIFE OF REAGENTS

Sealed / Unopened Reagents	
All sealed / unopened kit components are stable at 2-8 °C until the expiration date printed on the labels.	
Opened / Reconstituted Reagents	
Microtiter Plate	Return unused strips immediately to the aluminium pouch containing the desiccant packs and reseal along the entire edge of zip-seal. Store for up to 4 months at 2-8 °C.
Diluted Wash Buffer	Store for up to 4 months at 2-8 °C.
Calibrator	Store for up to 1 month at 2-8 °C. Do not freeze!
Controls	
Incubation Buffer	Store at 2-8 °C until expiration date printed on the labels.
Enzyme Labels	
TMB Substrate	
Stop Solution	Store at 18-28 °C until expiration date printed on the labels.

Table 2

PRECAUTIONS

Safety precautions

- Both, calibrator (B-GCO-CA) and controls (B-GCO-CONSET) of this kit contain components of human origin. Although tested and found negative for HBV surface antigen, HCV and HIV1/2 antibodies, the reagents should be handled as if capable of transmitting infections and should be handled in accordance with good laboratory practices using appropriate precautions.
- **Stop solution:** The stop solution (B-STTS) contains sulfuric acid (0.25 M). The reagent is an irritant to eyes, skin and mucous membranes. Avoid contact with eyes, skin and clothes. After contact with eyes or skin, wash immediately with plenty of water.
- **Reagents:** Avoid contact of reagents with the skin, eyes, or mucous membranes. If contact does occur, immediately wash with generous amounts of water; otherwise, irritation / burns can occur.
- Unused solution should be disposed of according to local state and federal regulations.

Technical precautions

- Read carefully the instructions prior to carrying out the test. Test performance will be adversely affected, if reagents are incorrectly diluted, modified or stored under conditions other than those as detailed in this instruction for use.
- Residues in the microtiter plate wells result from the production process. They are removed in the washing step (assay procedure step 3) and do not affect the results.
- Prepare reagents before starting the assay procedure. Reagents used in steps 3-9 must be cold (2-8 °C) and kept cold while pipetting and washing. Put the TMB substrate at room temperature (18-28 °C).
- Steps 3-9: Use cold (2-8 °C) reagents for all these steps and keep them cold while pipetting. Recommendation: Prepare the wash buffer the evening before performing the assay and place it into the fridge overnight.
- Wash steps 3, 6 and 9: The wash steps are crucial for removing residues in the microtiter plate wells resulting from the production process (step 3) as well as any unbound antibodies (steps 6 and 9).
 - Always perform the wash steps with cold (2-8 °C) wash buffer.
 - Make sure that all wells are completely empty after the last washing cycle.
- Step 9: Adjust TMB substrate to room temperature (18-28 °C) before using it.
- Step 11: Shake the microtiter plates during the incubation with substrate. Depending on the orbital plate shaker, we recommend 400-600 rpm. The solution should move in the wells but must not spill over.
- If an automated washer is used, "plate mode" should be chosen so that dispensing is performed sequentially on all strips before aspirating.
- Components must not be used after the expiry date printed on the labels.
- Do not mix different lots of reagents.
- Every effort should be made to ensure that no cross contamination occurs between reagents, samples or between wells.
- Microwells cannot be re-used.

MATERIALS REQUIRED BUT NOT PROVIDED

- Precision pipettes with disposable tips: 20 µL, 100 µL and 1000 µL pipettes.
- Disposable polystyrene or polypropylene tubes for the preparation of sample dilutions.
- 1000 mL cylinder for the reconstitution of the wash buffer.
- Squeeze bottle for wash buffer or automatic microtiter plate washer.
- Blotting paper.
- Orbital shaker for microtiter plates.
- Microtiter plate reader for measurement of absorbance at 450 nm.

SPECIMEN COLLECTION AND STORAGE

- The procedure requires <0.1 mL of blood or <50 µL of serum, respectively.
- Refer to page 6 to learn about the interference of haemolyzed, lipemic or icteric samples.
- Collect blood into plain tubes (no anti-coagulant), avoid haemolysis, leave to clot for one hour, centrifuge for 10 minutes at approximately 1500 x g at room temperature (18-28 °C), collect the serum.
- We recommend freezing aliquots of samples if you need to store samples in order to avoid repeated freezing / thawing.
- Store serum samples at ≤-20 °C up to 4 months. For long-term storage we recommend -70 °C (samples are stable for >1 year). Frozen samples should be thawed and vortexed thoroughly prior to use.

ASSAY PROCEDURE

You can choose between three basic options:

- (1) Detection of IgG/IgM mix-isotypes: step 4a-4e and 7
- (2) Detection of IgG and IgM isotypes: step 4a'-4f' and 7'
- (3) Two step approach: Option 1 → autoimmune antibody positive samples → Option 2.

Note: Equilibrate TMB substrate to room temperature (18-28 °C)

1. Dilute all samples to be investigated 1:50 with incubation buffer. Use 30 µL of serum + 1470 µL (cold: 2-8 °C!) incubation buffer. Mix by vortexing and leave diluted samples as well as reconstituted calibrator and controls for 30 minutes at 2-8 °C prior to pipetting (refer to step 4a and 4b).
2. Prepare a plate-frame with the required number of strips to test the samples. Reseal the remaining strips in the foil pouch together with the desiccant packs immediately. Store refrigerated.
3. Wash coated wells twice using at least 300 µL of cold! wash buffer per well. Empty wells and tap plate firmly onto blotting paper to remove remaining liquid completely.

Note: Use cold reagents in steps 3 to 9.

Option 1: Detection of IgG/IgM mix-isotypes

- 4a. Calibrator: Pipet 100 µL of calibrator into the well A1 (refer to figure 1A).
 - 4b. Controls: Pipet 100 µL of medium control into well B1, of low control into well A2 and of negative control into well B2 (refer to figure 1A).
- Note: If more than three strips per run are used, calibrator and controls can be tested in duplicates (see figure A).*
- 4c. Serum: Pipet 100 µL of diluted serum 1 into wells C1-E1 (refer to figure 1A).
 - 4d. Serum: Pipet 100 µL of diluted serum 2 into wells F1-H1 (refer to figure 1A).
 - 4e. Pipet 100 µL of diluted sera 3-24 into subsequent wells (refer to figure 1A).

Option 2: Detection of IgG isotypes

4a'. Calibrator: Pipet 100 µL of calibrator into the well A1 (refer to figure 1B).

4b'. Controls: Pipet 100 µL of medium control into well B1, of low control into well A2 and of negative control into well B2 (refer to figure 1B).

Note: If more than three strips per isotype are used, calibrator and controls can be tested in duplicates (see figure 1B).

4c'. Serum: Pipet 100 µL of diluted serum 1 into wells C1-E1 (refer to figure 1B).

4d'. Serum: Pipet 100 µL of diluted serum 2 into wells C1-H1 (refer to figure 1B).

4e'. Pipet 100 µL of diluted sera 3-12 into subsequent wells.

Detection of IgM isotypes.

4f'. Repeat steps 4a'-4e' using subsequent wells or a new microtiter plate if necessary (refer to figure 1B).

For options 1 and 2: Sample incubation and washes

5. Cover the plate with a plate sealer and incubate for 2 hours ±5 minutes at 2-8 °C (do not shake the plate).

6. Remove plate sealer. Empty the wells and wash three times using at least 300 µL of cold wash buffer (2-8 °C) per well. Empty the wells and strike the plate firmly onto blotting paper in order to remove washing buffer completely.

For option 1: Detection of IgG/IgM mix-isotype

7. Add 100 µL of enzyme label IgG/IgM mix to the wells.

For option 2: Detection of IgG and IgM isotypes

7'. Add 100 µL of enzyme label IgG or IgM to the respective wells (refer to figure 1B).

For option 1 and 2: Incubation with enzyme labels, washes, detection

8. Cover the plate with a plate sealer and incubate for 2 hours ±5 minutes at 2-8 °C (do not shake the plate).

9. Remove plate sealer. Empty the wells and wash three times using at least 300 µL of cold wash buffer (2-8 °C) per well. Empty the wells and strike the plate firmly onto blotting paper.

Note: Adjust TMB substrate solution to room temperature (18-28 °C).

10. Add 100 µL of TMB substrate solution to each well.

11. Cover plate with a plate sealer, incubate plate on an orbital plate shaker at 400-600 rpm for 30 ±2 minutes at 18-28 °C. Protect the plate from direct light.

12. Add 100 µL of stop solution to all wells. Proceed to step 13 within 30 minutes.

13. Read absorbance at 450 nm in a microtiter plate reader.

STANDARDIZATION

The calibrator included in this kit has been calibrated against internal reference material. It has been adjusted to 100 % ratio.

RESULTS AND CALCULATION

Calculation of results:

1. Record absorbance (OD) at 450 nm for each well (calibrator, controls and samples).
2. Average the duplicate calibrator and control values (if available).
3. Results are expressed as ratio of absorbance of samples and the (averaged) absorbance of the calibrator.

IgG/IgM Mix isotypes

$$\% \text{ Ratio: } \frac{\text{absorbance of samples or controls}}{\text{absorbance of calibrator}} \times 200$$

IgG and IgM isotypes

$$\% \text{ Ratio: } \frac{\text{absorbance of samples or controls}}{\text{absorbance of calibrator}} \times 100$$

Programs to calculate results as % ratio are available on most microplate readers.

Note: Results presented in tables 5 and 6 are examples. Calibrator and controls must be used in each individual assay.

QUALITY CONTROL

A good understanding of this instruction for use is necessary to obtain reliable results. These will be obtained only by using precise laboratory techniques (current GLP guidelines) and accurately following the instruction for use. Since there is no control serum for anti-ganglioside antibodies commercially available, we recommend using a positive, and negative serum pool for internal quality control. The calibrator must be within the established OD range. All controls must be within established confidence ranges (% ratio). The confidence ranges of the calibrator and for the controls are lot-specific and printed on the QC data sheet delivered with this kit.

Performance characteristics should be within established limits. If these characteristics are not in conformity with established limits and repetition excludes handling failures, check the following issues: i) Have all reagents, used in step 3-10, been kept at 2-8 °C? ii) accuracy of pipets, thermometers, and timers, iii) settings of ELISA washer and reader, iv) expiration date of the reagents v) storage and incubation conditions vi) colour of the TMB substrate solution (should be colourless) and vii) purity of the water.

PERFORMANCE LIMITATIONS

- The BÜHLMANN GanglioCombi™ *Light* ELISA has not been validated for plasmapheresis samples.

PERFORMANCE CHARACTERISTICS

Intra-assay precision (within-run): 1.4-5.7 %CV (mean: 3.4 %). For each of the three gangliosides coated on the microtiter plate of the BÜHLMANN GanglioCombi™ *Light* ELISA test two anti-ganglioside positive sera were selected. The serum samples were assayed in twelve replicates in a single run for each of the enzyme labels: IgG/IgM Mix, IgG and IgM. The intra-assay precision was determined to be in the range of 2.6 to 4.5 % CV for the enzyme label Mix, 2.5 to 5.7 % CV for the enzyme label IgG and 1.4 to 5.2 % CV for the enzyme label IgM. The study was performed with one reagent lot. The results are summarized in table 7, 8 and 9.

Inter-assay precision (between-run): 8.0-21.0 %CV (mean: 12.1 %). For each of the three gangliosides coated on the microtiter plate of the BÜHLMANN GanglioCombi™ *Light* ELISA test two anti-ganglioside positive sera were selected. The serum samples were assayed in single replicates in 20 independent runs, with one run per day, for each of the enzyme labels: IgG/IgM Mix, IgG and IgM. The inter-assay precision was determined to be in the range of 8.2 to 15.2 % CV for the enzyme label Mix, 11.6 to 21.0 % CV for the enzyme label IgG and 8.0 to 15.6 % CV for the enzyme label IgM. The study was performed with one reagent lot. The results are summarized in table 10, 11 and 12.

Limit of Blank (LoB): ≤6.1 %. The LoB was established according to the CLSI guideline EP17-A. Twelve blank replicates (incubation buffer) per ganglioside were assayed in a single run for all three enzyme labels: IgG/IgM Mix, IgG and IgM. The Limit of Blank (LoB) expressed as the % ratio to the calibrator absorbance was calculated to be 6.1 % for the enzyme label IgG/ IgM Mix, ≤3.1 % for the enzyme label IgG and ≤2.6 % for the enzyme label IgM. The highest LoB obtained with the three different enzyme labels was taken to determine the overall Limit of Blank (LoB). The study was performed with one reagent lot. The LoB was evaluated by using parametric analysis.

Limit of Detection (LoD): ≤8.1 %. The LoD was established according to the CLSI guideline EP17-A. For each of the three gangliosides coated on the microtiter plate a single clinical sample representing low antibody concentration was measured in twelve replicates, in a single run for each of the three enzyme labels: IgG/IgM Mix, IgG and IgM. The LoD expressed as the % ratio to the calibrator absorbance was calculated to be ≤8.1 % for the enzyme label IgG/ IgM Mix, ≤6.3 % for the enzyme label IgG, and ≤3.0 % for the enzyme label IgM. The highest LoD obtained with the three different enzyme labels was taken to determine the overall Limit of Detection (LoD). The study was performed with one reagent lot. The LoD was calculated using parametric analysis.

Functional sensitivity: ≤ 8.1 %. A precision profile was generated from the results obtained in the inter-assay precision study described above, for two positive sera samples for each of the three gangliosides coated on the microtiter plate and detected with each of the enzyme labels: IgG/IgM Mix, IgG and IgM. The data points were fitted using a polynomial non-linear fit with a cubic function.

The 95 %-confidence intervals of the fit were determined and the functional sensitivity was determined as the intersection of the cubic fit with the 20 % CV acceptance criterion. The results are summarized in figure 2.

Linearity: The linear range of the BÜHLMANN GanglioCombi™ *Light* ELISA was determined according to the CLSI guideline EP06-A. Multiple sera with concentrations over the entire measuring range of the test, allowing the evaluation of the majority of ganglioside / enzyme label combinations, were used. Serum samples were diluted according to the instruction for use. In case of high positive serum samples a higher dilution of 1:2000 was applied. Subsequently dilution series from each sample were prepared in graduations of 10 % using incubation buffer as the diluent. Linearity was defined as the interval in which the relative difference between the linear and, if significant, higher order polynomial fit was below 20 %. For ratios ≤25 % an absolute difference of below 5 % ratio was allowed. A linear range covering and exceeding the BÜHLMANN GanglioCombi™ *Light* ELISA titre categories of <30 % - negative; 30-50 % - grey-zone; 50-100 % positive and >100 % strongly positive was confirmed. The results are summarized in table 13.

INTERFERING SUBSTANCES

No interference is detected with the following substances up to the following concentrations: unconjugated bilirubin (icteric sera): 40 mg/d; conjugated bilirubin (icteric sera): 60 mg/dL; haemoglobin (haemolyzed sera): 400 mg/dL; haemolyzate (haemolyzed blood): 400 mg/dL and Triglycerides (Intralipid®): 3000 mg/dL.

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