# Extraction Kit

for BÜHLMANN fCAL<sup>®</sup> Assay

B-CAL-EX12

Revision date: 2015-11-05

## ENGLISH

### INTENDED USE

The extraction buffer (B-CAL-EX) is a reagent designed for the extraction of calprotectin in faeces using the BÜHLMANN Smart-Prep or Manual weighing method.

### REAGENTS SUPPLIED

Reagents	Quantity	Code	Preparation
Extraction Kit extraction buffer	12 bottles x 125 mL	B-CAL-EX12	Ready to use

Table 1

### MATERIALS REQUIRED BUT NOT PROVIDED

- Smart-Prep (Code : B-CAL-RD)
- 10 µL disposable inoculation loops
- 15 mL polypropylene tubes with screw caps required for standard extraction procedure
- Multi tube vortex mixer
- Precision balance (10-150 mg)
- Micro centrifuge (≥3000 g)

### REAGENT STORAGE AND STABILITY

- **Do not freeze reagents!** Shelf life of unopened/opened reagents at 2- 8 °C: see expiry date on package label.

### WARNINGS AND PRECAUTIONS

- Read carefully the instructions prior to use the buffer. Buffer performance will be adversely affected, if the reagent is incorrectly stored as detailed in this instruction for use.
- Components must not be used after the expiry date printed on the labels.
- Let the reagent adjust to reach room temperature. Mix well (vortex) the reagent before use.

### Extraction

- To receive quantitative results it is important to homogenize the stool sample completely with the extraction buffer.

### SPECIMEN COLLECTION AND HANDLING

If the extraction devices are used less than 1 g of native stool sample is needed for the extraction procedure. Collect stool samples into plain tubes. The samples can be stored refrigerated at 2-8°C for at least 6 days. The extracts are stable for at least 7 days at 2-8°C and for at least 24 months at -20°C.

**Important:** The sample must be collected without any chemical or biological additions in the collection device.

### STOOL SAMPLE EXTRACTION

#### with BÜHLMANN Smart-Prep

How to use the Smart-Prep (Code: B-CAL-RD) see the corresponding Instruction for Use (IFU) of the device.

#### with original<sup>1)</sup> Manual weighing method

1. Label and weigh (tare) the empty polypropylene tube together with the inoculation loop.
2. Take out 50 to 100 mg of the stool sample by means of the inoculation loop and place it into the pre-weighted tube.
3. Calculate the net amount of sample, break off the inoculation loop and leave the lower part of the loop in the tube.

4. Add Extraction Buffer according to the formula  
 $x \text{ mg stool} \times 49 = y \text{ } \mu\text{L extraction buffer}$   
(e.g.  $50 \text{ mg stool} \times 49 = 2450 \text{ } \mu\text{L buffer}$ )  
to the tube and close the tube
5. Homogenize the sample on a multi tube vortexer by vigorous shaking (at highest speed) for 30 minutes.
6. Transfer the homogenate into a 2 mL Eppendorf tube and centrifuge in a microcentrifuge for 5 minutes at 3'000 x g.
7. Take the supernatant into a fresh, labelled tube and continue with the ASSAY procedure of the corresponding BÜHLMANN kit.

**Important:** The extraction with the buffer (B-CAL-EX) is only useful in combination with the BÜHLMANN Kit

### SHIPPING DAMAGE

Please notify your distributor, if this product was received damaged.

### REACH

None of the materials and reagents in the kit require a Safety Data Sheet (MSDS) according to CLP-Regulation (EC) No 1272/2008 and directive EC 1907/2006 (REACH).

<sup>1)</sup> according BÜHLMANN fCAL<sup>®</sup> ELISA test

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