



BÜHLMANN GanglioCombi® MAG ELISA

Pipetting Template



Pipetting Guidance

Precoated Microtiter Plate



wash 2x

100 µL of Calibrator, Controls or Serum Samples (1:50)



incubate 2 hours ± 5 min at 2-8°C wash 3x
(use cold wash buffer from refrigerator)

add 100 µL of Enzyme Label(s)



incubate 2 hours ± 5 min at 2-8°C wash 3x
(use cold wash buffer from refrigerator)

add 100 µL of TMB Substrate
(ambient temperature)


incubate 30 min ± 2 min at 18-28°C



add 100 µL of Stop Solution

Read absorbance at 450 nm (within 30 minutes)

Autoimmune Peripheral Neuropathy Related Research: Three Easy Steps for Complete Assessment



*„MAG“ is a chemically synthesized, sulfated disaccharide which mimics a MAG (HNK-1) epitope.



BÜHLMANN Laboratories AG
Switzerland
info@buhlmannlabs.ch
www.buhlmannlabs.ch
Phone +41 61 487 12 12

BÜHLMANN Diagnostics Corp
USA / Canada
info@buhlmannlabs.com
www.buhlmannlabs.com
Phone 844 300 9799

BÜHLMANN GanglioCombi® MAG ELISA (EK-GCM), Anti-MAG Autoantibodies ELISA (EK-MAG), anti-GM1 Autoantibodies ELISA (EK-GM1-GM) and anti-SGPG Autoantibodies ELISA (EK-SGPG) are Health Canada Licensed.

BÜHLMANN GanglioCombi® MAG ELISA anti-MAG Autoantibodies ELISA κ/λ Reagent Set anti-SGPG Autoantibodies ELISA BÜHLMANN GanglioCombi® Light ELISA anti-GM1 Autoantibodies ELISA

EK-GCM-U
EK-MAG-U
B-KLSET-U
EK-SGPG-U
EK-GCL-S-U
EK-GM1-GM-U

BÜHLMANN GanglioCombi® is a registered trademark of BÜHLMANN in Switzerland.