anti-MAG

anti-Myelin Associated Glycoprotein Autoantibodies ELISA

EK-MAG 96 tests

Revision date: 2012-11-16

ENGLISH

INTENDED USE

The BÜHLMANN anti-MAG ELISA is intended for the quantitative *in vitro* diagnostic determination of human IgM-antibodies directed against Myelin Associated Glycoprotein (MAG) (1, 2).

PRINCIPLE OF THE ASSAY

The anti-MAG autoantibodies ELISA employs the quantitative enzymatically amplified sandwich-type immunoassay technique. HIGHLY PURIFIED MAG FROM HUMAN BRAIN (2) has been precoated onto a microtiter plate. Calibrators and patient sera are incubated for two hours in the microtiter wells and any anti-MAG autoantibodies present are bound by the immobilized human MAG. After washing away any unbound substances, horseradish peroxidase (HRP) labeled antibodies against human IgM are added to the wells and incubated for another two hours. After a wash step, the substrate solution containing tetramethylbenzidine (TMB) is added to the wells and incubated for 30 minutes. A blue coloration develops in proportion to the amount of anti-MAG autoantibodies bound in the initial step. The color development is stopped by adding an acidic stop solution (H₂SO₄) which turns the blue solution to yellow. The intensity of the color absorbance is measured in a microtiter plate reader at a wavelength of 450 nm. The absorbance measured is directly proportional to the concentration of anti-human MAG autoantibodies. A set of human anti-MAG autoantibody calibrators is used to plot a standard curve of absorbance versus human anti-MAG autoantibody titer units from which the concentrations of anti-human MAG autoantibodies in the unknowns can be calculated.

REAGENTS SUPPLIED AND PREPARATION

Reagents	Quantity	Code	Reconstitution
Microtiter Plate 96 wells precoated with human MAG	12 x 8-well strips with holder	B-MAG-MP	Ready to use
Plate Sealer	3 pieces		
Wash Buffer Concentrate (10x)	1 bottle 100 ml	B-MAG-WB	Dilute with 900 ml of deionized water
Incubation Buffer with preservatives	1 bottle 100 ml	B-MAG-IB	Ready to use
Calibrators A to D¹ Human serum with preservatives	4 vials	B-MAG-CASET	Add 1 ml of Incubation Buffer
Low and High Control² Human serum with preservatives	2 vials	B-MAG-CONSET	Add 1 ml of Incubation Buffer
Enzyme Label IgM Anti-human IgM antibody conjugated to HRP in a protein-based buffer with preservatives	1 vial 11 ml	B-MAG-ELM	Ready to use Blue solution
TMB Substrate TMB in Citrate buffer with Hydrogen Peroxide	1 vial 11 ml	B-TMB	Ready to use
Stop Solution 0.25 M Sulfuric acid	1 vial 11 ml	B- STS	Ready to use Corrosive agent

Table 1

¹ After reconstitution, Standards A, B, C and D contain 70000, 15000, 3000 and 1000 Böhmann Titer Units (BTU) of anti-MAG antibodies, respectively.

² Controls contain lot-specific amounts of anti-MAG antibodies. Refer to the additional control data sheet for exact concentrations.

STORAGE AND SHELF LIFE OF REAGENTS

Unopened Reagents	
Store at 2-8°C. Do not use past kit expiration date printed on the label.	
Opened / Reconstituted Reagents	
Microtiter Plate	Return unused strips immediately to the foil pouch containing the desiccant packs and reseal along the entire edge of zip-seal. Store for up to 2 months at 2-8°C.
Wash Buffer diluted	Store for up to 2 months at 2-8°C.
Calibrators	Store for up to 2 months at -20°C.
Controls	
Incubation Buffer	Store at 2-8°C until expiration date.
Enzyme Label IgM	
TMB Substrate (protect from light)	
Stop Solution	Store at 18-28°C until expiration date.

Table 2

PRECAUTIONS

SAFETY PRECAUTIONS

- Both, Calibrator (B-MAG-CASET) and Controls (B-MAG-CONSET) of this kit contain components of human origin. Although tested and found negative for HBV surface antigen, HCV and HIV1/2 antibodies, the reagents should be handled as if capable of transmitting infections and should be handled in accordance with good laboratory practices using appropriate precautions.
- Substrate and Stop Solution:** The Substrate Solution (B-TMB) contains Tetramethylbenzidine (TMB), hydrogen peroxide and dimethylformamide. The Stop Solution (B-STs) contains sulfuric acid. Each of those reagents is irritant to eyes, skin and mucous membranes. Avoid contact with eyes, skin and clothing. After contact with eyes or skin, wash immediately with plenty of water.
- Unused solution should be disposed of according to local State and Federal regulations.

TECHNICAL PRECAUTIONS

Kit components

- Read carefully the instructions prior to carrying out the test. Test performance will be adversely affected, if reagents are incorrectly diluted, modified or stored under conditions other than those as detailed in this instruction for use:
- Residues in the microtiter plate wells** result from the production process. They are removed in the washing step (Assay procedure step 3) and do not affect the results.
- Steps 3-9:** Use cold (2-8°C) reagents for all these steps and keep them cold while pipetting.
- Steps 3, 6, 9:** Make sure that the wells are completely empty after the last washing cycle.
- Step 9:** Adjust TMB Substrate to room temperature (18-28°C) before using it.
- Step 11:** Shake microtiter plates during the incubation with substrate. Depending on the plate shaker, we recommend 400-600 rpm. The solution should be moved in the wells but must not spill over.
- If an **automated washer is used**, "plate mode" should be chosen so that dispensing is performed sequentially on all strips before aspirating.
- Components must not be used after the expiry date printed on the labels.
- Do not mix different lots of reagents.

- Every effort should be made to ensure that no cross contamination occurs between reagents, samples or between wells.
- Microwells cannot be re-used.

MATERIALS REQUIRED BUT NOT PROVIDED

- Precision pipettes with disposable tips: 2 µl, 100 µl and 1 ml pipettes.
- Disposable polystyrene or polypropylene tubes for the preparation of sample dilutions.
- 1000 ml cylinder for the dilution of the Wash Buffer Concentrate.
- Microtiter plate washer or squeeze bottle for Wash Buffer.
- Microtiter plate rotator.
- Microtiter plate reader for measurement of absorbance at 450 nm.

SPECIMEN COLLECTION AND STORAGE

The procedure calls for <0.1 ml of blood or <50 µl of serum. Lipemic, hemolytic and icteric samples should not be used in this assay. Lipemic samples can be avoided by asking patients to fast for at least 12 hours prior to the sample being taken. Collect blood into plain tubes, avoid hemolysis, leave to clot for one hour at RT (18-28°C), centrifuge for 15 minutes at approximately 1800 x g at RT and collect the serum.

Store serum samples at ≤-20°C. Samples are stable for ≥1 year if stored at ≤-20°C. Avoid repeated freeze-thaw cycles. Frozen samples should be thawed and mixed thoroughly by gentle swirling or inversion prior to use.

ASSAY PROCEDURE

1. Dilute all patient samples 1:1000 with Incubation Buffer (e.g. 2 µl of serum + 2 ml of Incubation Buffer). Allow diluted samples to set for one hour at 18-28°C vortex from time to time. Put samples for 10 minutes on ice prior to pipetting in step 4c.
 2. Prepare a plate with sufficient strips to test the desired number of calibrators, controls and samples. Remove excess strips from the holder and reseal them in the foil pouch together with the desiccant packs **without delay**. Store refrigerated.
- Important: Use refrigerated reagent solutions in steps 3. to 9.**
3. Wash the coated wells four times using at least 300 µl of refrigerated Wash Buffer per well. Empty the wells and strike the plate firmly onto blotting paper.
 - 4a. Pipet 100 µl of Incubation Buffer (as blank) in duplicate into wells A1+A2.
Pipet 100 µl of Calibrator A in duplicate into wells B1+B2.
Pipet 100 µl of Calibrator B in duplicate into wells C1+C2.
Pipet 100 µl of Calibrator C in duplicate into wells D1+D2.
Pipet 100 µl of Calibrator D in duplicate into wells E1+E2.
 - 4b. Pipet 100 µl of the Low Control in duplicate into wells F1+F2.
Pipet 100 µl of the High Control in duplicate into wells G1+G2.
 - 4c. Pipet 100 µl of each diluted sample in duplicate into the subsequent wells.
 5. Cover the plate with a Plate Sealer and incubate for 2 hours (± 5 min) at 2-8°C.
 6. Remove and discard the Plate Sealer. Empty the wells and wash four times using at least 300 µl of

refrigerated Wash Buffer per well. Empty the wells and strike the plate firmly onto blotting paper.

7. Add 100 µl of Enzyme Label IgM to all wells.
8. Cover the plate with a Plate Sealer and incubate for 2 hours (± 5 min) at 2-8°C.
9. Remove and discard the Plate Sealer. Empty the wells and wash four times using at least 300 µl of **refrigerated** Wash Buffer per well. Empty the wells and strike the plate firmly onto blotting paper.

Important: Allow TMB substrate solution to reach 18-28°C.

10. Add 100 µl of the TMB Substrate Solution to each well.
11. Cover the plate with a Plate Sealer, place the plate on a plate mixer set at 800-1000 rpm, protect the plate from direct light and incubate for 30 minutes (± 5 min) at 18-28°C.
12. Add 100 µl of Stop Solution to all wells. Remove air bubbles with a pipette tip. Proceed to step 13. within 30 minutes.
13. Read the absorbance at 450 nm in a microtiter plate reader. If wavelength correction is available, set instrument to dual wavelength measurement at 450 nm with background wavelength correction set at 600 or 620 nm.

STANDARDIZATION

The Calibrators included in this kit has been calibrated against internal reference pool. The reference pool consists of more than 10 human sera containing low, medium and high titers of anti-MAG antibodies, respectively. Titers and specificity of each serum in the reference pool were first analyzed by anti-MAG immunoblots.

The **Bühlmann Titer Units (BTU)** were established as follows:

- Normal donor samples were assayed according to the anti-MAG ELISA assay procedure.
 - Serially diluted samples of the Reference Pool were assayed in the same run.
 - The dilution at which the Reference Pool sample falls short of the cut-off corresponds to the titer of the reference pool.
- NOTE: Serum titer values depend on the assay method and, in particular, on the specificity and the cut-off values established with an individual assay method. Therefore, titer values obtained with different methods cannot be compared directly.

RESULTS AND CALCULATION

Standard Curve

- Record the absorbance at 450 nm for each calibrator.
- Average the duplicate values.
- Plot the average absorbance values (vertical y-axis) versus BTU (see next chapter) of the calibrators (horizontal x-axis) using a lin/log graph paper.
- Draw the best fitting curve or calculate the standard curve using a spline smoothed fitting algorithm.

Samples and Controls

- Record the absorbance at 450 nm for each sample and control well.
- Average the duplicate values.
- Locate the absorbance values of the samples and controls on the vertical axis, draw a horizontal line intersecting the standard curve and read from the horizontal axis.

NOTE: If the microtiter plate reader is not capable of reading absorbance greater than 2 or greater than the absorbance of the highest calibrator (Calibrator A), a second reading at a wavelength of 490 or 492 nm is recommended (reference filter at 600 or 620 nm if available). In this case, proceed to construct a second standard curve with the absorbance readings of all calibrators at 490 or 492 nm. The

concentration of the off-scale samples at 450 nm are then read from the new standard curve as described above. The readings at 490 or 492 nm should not replace the on-scale readings at 450 nm.

Note: The results presented in Table 11 and Figure 1 are examples. Calibrator and Controls must be used in each individual assay.

These results and standard curve are provided for demonstration purposes only. A standard curve must be generated for each set of samples to be assayed.

QUALITY CONTROL

A thorough understanding of this package insert is necessary for the successful use of the product. Reliable results will be obtained only by using precise laboratory techniques (current GLP guidelines) and accurately following this package insert.

Since there is no control serum for anti-MAG antibodies commercially available, we recommend to use a positive serum pool for internal quality controls.

All controls must fall within established confidence limits. The confidence limits for the Controls are lot-specific and printed on the additional data sheet.

The reproducibility of standard curve parameters and control values should be within established limits of laboratory acceptability. If the precision of the assay does not correlate with the established limits and repetition excludes errors in technique, check the following issues: i) pipetting, temperature controlling and timing devices ii) ELISA reader settings iii) expiration dates of reagents iv) storage and incubation conditions v) TMB Substrate Solution should be colorless vi) purity of water.

PERFORMANCE LIMITATIONS

Test results should be interpreted in conjunction with information available from clinical assessment of the patient and other diagnostic procedures.

PERFORMANCE CHARACTERISTICS

Intra-Assay Precision (Within-Run): 6.5%. The intra-assay precision was calculated from the results of 20 pairs of values obtained in a single run (*cf.* Table 13).

Inter-Assay Precision (Run-to-Run): 15.4%. The inter-assay precision was calculated from the results of 20 pairs of values obtained in 20 different runs (*cf.* Table 14).

Dilution Linearity/Parallelism: 147%. 7 human serum samples containing high titer of anti-MAG antibodies were diluted with Incubation Buffer 1:1000 to 1:64000, left for one hour at 18-28°C and subsequently assayed according to the assay procedure (*cf.* Table 15). It is suggested that the relatively high deviation in about 50% of the samples is due to antibody aggregations. In general, pathological sera show strongly elevated titer of autoantibodies therefore it has no influence to the positive/negative discrimination.

Limit of Blank (LoB): 444 BTU. 20 duplicates of incubation buffer were assayed in a single run. Mean and standard deviation were calculated for the absorbance values. The minimal detectable dose of anti-MAG antibodies was calculated to be 444 BTU by adding two standard deviations to the mean absorbance and intersecting this value with the standard curve obtained in this run.

Limit of Quantification (LoQ): 900 BTU. 20 duplicates of low titer anti-MAG autoantibodies were assayed in a single run. Mean, standard deviation (SD) and coefficient of

variation (CV) were calculated from the absorbance values. We found a titer of 900 BTU with a CV less than 10%.

Specificity: Four sets of experiments were performed to assess the specificity of the BÜHLMANN anti-MAG ELISA:

1. **NEUTRALIZATION OF ANTI-MAG AUTOANTIBODIES:** Five sera with high anti-MAG titers could be inhibited from binding to the microtiter plates coated with MAG in a concentration-dependent manner when preincubated for one hour with Incubation Buffer supplemented with 1 to 200µg/ml of MAG prior to testing in the ELISA.

2. **SPECIFICITY OF ANTI-MAG AUTOANTIBODY BINDING:** Five sera with medium to high anti-MAG autoantibody titers and five negative sera were tested on a GanglioCombi (EK-GCO) plate coated with asialo-GM1, GM1, GM2, GD1a, GD1b and GQ1b. None of the sample showed ratio higher than 10% (Cut-off ratio >12%). Similarly, six sera with high anti-Ganglioside autoantibody titers (see above) and five negative sera, were tested on plates coated with MAG. None of these disease state sera gave a signal higher than 300 BTU.

3. **COMPARISON TO IMMUNOBLOT (WESTERN BLOT):** 127 patient sera (40 female; 87 male; age 4-90 years) with neurological diseases suspectedly caused by anti-MAG autoantibodies were tested by the ELISA and a western blot method, respectively. MAG isolated from human central nervous system (CNS) was used in both methods. All 12 sera showing medium to very high titers of anti-MAG IgM autoantibodies in the ELISA procedure were also clearly positive when analyzed on immunoblots. 114 out of 115 sera showing anti-MAG IgM antibody titers below the ELISA cut-off value were also negative when analyzed on immunoblots, whereas 1 serum was slightly positive by the immunoblot method only (F. Ferracin and A.J. Steck, unpublished results).

4. **COMPARISON TO INDIRECT IMMUNOFLUORESCENCE ASSAY (IFA):** Sera from 150 patients (unknown sex and age) diagnosed as having IgM monoclonal gammopathies were tested by the ELISA and an IFA method, respectively. In the IFA method, sera were incubated on frozen and acetone-fixed monkey sciatic nerve sections and bound anti-MAG antibodies were detected by FITC-labeled anti-human IgM antibody. 26 patients (17.3%) were positive and 115 patients (76.7%) were negative in both methods, whereas 5 sera (3.3%) were positive in the IFA method only and 4 sera (2.7%) were only positive in the ELISA method, respectively (3).

REFERENCE INTERVALS AND CUT-OFF

The frequency of anti-MAG antibodies in normal human sera was determined using blood samples from asymptomatic volunteer blood donors (adult men and women at the age of 18 to 70 years). 150 samples were assayed according to the assay procedure and the results listed in Table 12 were obtained.

NOTE: These titer ranges should be used as guidelines only. It is recommended that each laboratory establishes its own expected ranges for its patient population.

Proposed Cut-off Titer

The mean + 3SD results in a technical cut-off value of 729 BTU. For practical reasons, we recommend to use a CUT-OFF VALUE OF 1000 BTU, corresponding to the lowest calibrator (= Calibrator D) in the standard curve.

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
Figure 1: **Example of Standard Curve**Table 12: **Cut-off values and Standardization**

Total (n)	150
Range (BTU)	0 -832
Mean Value (BTU)	173
SD (BTU)	186
Mean + 3 SD (BTU)	729

Table 13: **Intra-Assay Precision (Within-Run)**

Sample Type	Mean [BTU]	SD [BTU]	CV [%]
Serum 1 (Low)	3764	176	4.7
Serum 2 (High)	35384	2920	8.3
Mean			6.5

Table 14: **Inter-Assay Precision (Run-to-Run)**

Sample Type	Mean [BTU]	SD [BTU]	CV [%]
Serum 3 (Low)	4208	714	17.0
Serum 4 (High)	17112	2362	13.8
Mean			15.4

Table 15: **Dilution Linearity**

Serum Samples	Range Min – Max	Mean Observed/Expected
5	78 – 96%	86%
6	99 – 119%	107%
7	105 – 151%	123%
8	162 – 229%	206%
9	135 – 231%	196%
10	112 – 292%	233%
11	59 – 98%	76%
Mean		147%

Table description: cf. “Results”, and “Cut-Off Values and Standardization”) and “Performance Characteristics”

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
1. Quarles, R.H and Weiss, M.D.: Autoantibodies associated with peripheral neuropathy. *Muscle Nerve* 22, 800-822 (1999).
2. Burger, D. and Steck, A.J.: Neuropathies associated with anti-myelin antibodies. In: Hohlfeld, R. (Ed.): *Immunology of neuromuscular disease* 24, 7-32 (1994).
3. Schmit, P. et al.: Evaluation of four different methods for the detection of anti-myelin autoantibodies. Poster presented at the MEDLAB97, 12th IFCC, Europ Congr of Clin Chem, Basel. Abstract C139 (pp 229-230)

Anti-MAG ELISA

Precoated Microtiter Plate

TIME TO RESULT: 4.5 HOURS

APPENDIX IV
SYMBOLS/SYMBOLE/ SYMBOLES/SIMBOLI/ SIMBOLOS

Symbol	Explanation
	Use By Verwendbar bis Utiliser jusqu'au Utilizzare entro Fecha de caducidad
REF	Catalogue number Bestellnummer Référence du catalogue Numero di catalogo Número de catálogo
LOT	Batch code Chargenbezeichnung Code du lot Codice del lotto Codigo de lote
IVD	<i>In Vitro</i> Diagnostic Medical Device <i>In Vitro</i> Diagnostikum Dispositif médical de diagnostic <i>in vitro</i> Dispositivo medico-diagnostico <i>in vitro</i> Producto sanitario para diagnóstico <i>in vitro</i>
	Contains sufficient for <n> tests Ausreichend für "n" Ansätze Contenu suffisant pour „n“ tests Contenuto sufficiente per „n“ saggi Contenido suficiente para <n> ensayos
	Consult Instructions for Use- Gebrauchsanweisung beachten Consulter le mode d'emploi Consultare le istruzioni per l'uso Consulte las instrucciones de uso
	Temperature limitation Zulässiger Temperaturbereich Limites de température Limiti di temperatura Limite de temperatura
BUF INC	Incubation Buffer Inkubations-Puffer Tampon d'incubation Tampone di incubazione Tampón de incubación

Symbol	Explanation
MP	Microtiter Plate Mikrotiter-Platte Microplaque Micropiastra Microplaca
BUF WASH 10X	Wash Buffer Concentrate (10x) Wasch-Puffer Konzentrat (10x) Tampon de lavage concentré (10x) Tampone di lavaggio concentrato (10x) Tampón de lavado concentrado (10x)
CAL A - CAL D	Calibrator A - D Kalibrator A - D Calibreur A - D Calibratore A - D Calibrador A - D
CONTROL L	Low Control Kontrolle tief Contrôle bas Controllo basso Control bajo
CONTROL H	High Control Kontrolle hoch Contrôle élevé Controllo alto Control alto
EL IgM	Enzyme Label IgM Enzymmarker IgM Marqueur enzymatique IgM Marcato enzimatico IgM marcador enzimático IgM
SUBS TMB	TMB Substrate TMB Substrat Substrat TMB Substrato di TMB Substrato de TMB
SOLN STOP	Stop Solution Stopp-Lösung Solution stop Soluzione stoppante Solución de parada

Reg. Nr. 22176

Printing Date
2013-02-06