ANGIOTENSIN II

RIA

RK-A22 100 tests

Revision date: 2012-11-20

ENGLISH

INTENDED USE

The BÜHLMANN Angiotensin II is a double antibody radioimmunoassay designed for the quantitative *in vitro* diagnostic determination of immunoreactive Angiotensin II in extracted EDTA plasma.

PRINCIPLE OF THE ASSAY

The BÜHLMANN Angiotensin II radioimmunoassay measures immunoreactive Angiotensin II (Ang II) by a double antibody radioimmunoassay using a modification of the method reported by Emanuel *et al.* (1). Extracted EDTA plasma samples, calibrators and controls are first preincubated for 16 hours with an anti-Ang II antibody. ¹²⁵I-Ang II is added and competes with Ang II present in samples, calibrators and controls for the same antibody binding sites in a second 6 hours incubation step. After this incubation, a solid-phase second antibody is added to the mixture. The antibody-bound fraction is precipitated and counted in a gamma-counter.

It is recommended to treat plasma samples with bestatin or an equivalent angiotensinase inhibitor during blood sampling in order to prevent degradation of Ang II (5-8; see also Specimen Collection and Storage).

REAGENTS SUPPLIED AND PREPARATION

| Reagents | Quantity | Code | Reconstitution |
|---|------------------|--------------|---|
| Extraction Columns (1 ml) phenylsilylsilica reversed phase extraction columns | 10 pieces | B-AEC | |
| Tris Buffer | 1 vial 100 ml | B-A22-TB | Ready to use |
| Antiserum rabbit anti-Ang II antibody | 1 vial 10 ml | B-A22-AS | Ready to use |
| Tracer ¹²⁵ I-Ang II | 1 vial 11 ml | B-A22-TR | Ready to use |
| Calibrators¹⁾ lyoph. synthetic Ang II | 5 vials | B-A22-CASET | Reconstitute each vial with 5 ml of Tris Buffer |
| Controls Normal / High²⁾ lyophilized Ang II in a buffer matrix | 2 vials | B-A22-CONSET | Reconstitute each vial with 5 ml of Tris Buffer |
| Second Antibody Solid phase bound anti-rabbit second antibody | 1 vial 11 ml | B-A22-AB2 | Ready to use |

Table 1

¹⁾ After reconstitution, the calibrator solutions will contain 2, 5, 20, 100 and 500 pg/ml of Ang II, respectively.

²⁾ See QC data sheet for the lot specific amount of Ang II

STORAGE AND SHELF LIFE OF REAGENTS

| Unopened Reagents | |
|--|--|
| Columns should be stored at 18-28°C. All other kit components are stable at 2-8°C until the expiration date printed on the labels. | |
| Opened / Reconstituted Reagents | |
| Extraction Columns | Used columns should be stored at 18-28°C and protected from dust and light. |
| Tris Buffer | Stable at 2-8°C until expiration date printed on the label. |
| Antiserum | |
| Tracer | |
| Calibrators | Stable for at least 2 months after reconstitution at -20°C. |
| Controls | |
| Second Antibody | Store refrigerated (Do not freeze!) Stable at 2-8°C until expiration date printed on the label. |

Table 2

PRECAUTIONS

SAFETY PRECAUTIONS

- Radioactive Material: This kit contains radioactive material which does not exceed 56 kBq of ¹²⁵Iodine.
- Receipt, acquisition, possession, use, and transfer are subject to the local regulations. Unused solutions and radioactive waste should be disposed of according to local State and Federal regulations.

TECHNICAL PRECAUTIONS

- Read carefully the instructions prior to carrying out the test. Test performance will be adversely affected if reagents are incorrectly diluted, modified or stored under conditions other than those as detailed in this instruction for use.
- Incorrect results for standard curve, controls or samples may be obtained, if the 2nd antibody was not mixed sufficiently before pipetting. Do not freeze the 2nd antibody!
- Components must not be used after the expiry date printed on the labels.
- Do not mix different lots of reagents.
- Let the reagents adjust to reach room temperature. Reconstitute the lyophilized reagents as indicated. Mix well (vortex) the reagents before use.
- Counting time should be selected in order to keep statistical counting error small: e.g., at 2000 cpm the counting error is at 5%; at 10000 cpm it is only 1%.
- If the initial concentration of an unknown sample reads above the highest calibrator, the sample should be further diluted with Tris buffer and tested again according to the assay procedure.

MATERIALS REQUIRED BUT NOT PROVIDED

- 100 µl, 500 µl, 1000 µl, and 5000 µl precision pipettes with disposable tips.
- Disposable polypropylene tubes (13 x 100 mm) for preparation of plasma extracts.
- Conical polystyrene tubes for the assay (e.g. Sarstedt 57.477).
- Extraction vacuum manifold for applying the extraction columns (optional).
- Methanol (HPLC grade recommended).
- Distilled or deionized water.
- Bestatin solution¹⁾ for sample collection (optional).
- Ice-bath
- Refrigerated centrifuge.
- Vortex mixer.
- Stir bar and magnetic stirrer.
- Aspiration device.
- Evaporator.
- Gamma-counter.

¹⁾The Bestatin reagent can be ordered separately (order code: B-A22-BST).

SPECIMEN COLLECTION AND STORAGE

The procedure calls for 1.2 ml of EDTA plasma for duplicate determinations. Appropriate sample collection is essential to ensure accurate results of the Ang II analysis. Hemolyzed, highly icteric or lipemic samples may adversely effect results.

Recommended Collection Procedure: It is important to prevent the *in vitro* generation of Ang I and Ang II from angiotensinogen and Ang I by renin and ACE, respectively, as well as the *in vitro* degradation of Ang II to Ang III by angiotensinases. The *in vitro* generation and degradation of Ang II can be minimized by following this recommended collection procedure:

Collection: Collect 5 ml of blood from **recumbent** patient into a pre-chilled EDTA venipuncture tube or plastic syringe. Immediately add 100 µl of ready to use Bestatin solution (see above) or equivalent angiotensinase inhibitor (*cf.* 5,7-9) per 5 ml EDTA-blood sample. Place sample immediately in an ice-bath and proceed with plasma separation within **1 hour**.

Separation: Centrifuge sample for 15 minutes at 1000 x g and 2-8°C, separate the plasma from the cells, aliquot and freeze the specimen immediately in polypropylene tubes at -20°C or less or extract plasma samples without delay.

Alternative Collection Procedure: Collect 5 ml of blood from **recumbent** patient into a pre-chilled EDTA venipuncture tube or a plastic syringe and immediately place the sample on ice. Centrifuge at 1000 x g and 2-8°C **immediately** after collection. Separate the plasma from the cells, aliquot and freeze the specimen at -20°C in polypropylene tube or extract plasma samples without delay.

COLUMN EXTRACTION OF PLASMA SAMPLES

The use of a reversed phase extraction method is essential for obtaining accurate Ang II results. The best recoveries are obtained with phenylsilylsilica extraction columns (*cf.* 2-5). The extraction columns provided with this kit can **each be utilized up to five times** if used according to the extraction procedures described in this protocol. Used columns should be stored at 18-28°C and protected from light and dust. In order to avoid clogging of columns, **filter or centrifuge samples containing particles** such as fibrin clots prior to extraction. This extraction method usually results in recoveries greater than 90 % (*cf.* 5,10,12). The following extraction procedure was **tested and validated for human EDTA plasma** samples. If other specimens are used, it is recommended to validate the extraction recovery using a ¹²⁵I-Angiotensin II-spiked specimen (e.g. 1 ml of specimen spiked with 50 µl of ¹²⁵I-Angiotensin II).

Extraction Procedure Using Vacuum Manifold

With the help of negative pressure the fluids will be passed through the column. The procedure is the same as described in the Extraction procedure using centrifugation (see below).

The following flow rates should be used:

Sample application and elution should be done with a flow rate of 2 ml/min

All the other fluids can pass the column with a flow rate of 5 ml/min.

Extraction Procedure using Centrifugation

COLUMN PREPARATION AND CONDITIONING

- Mark one extraction column for each sample to be extracted and place into polypropylene centrifugation tubes.
- Add 2x 1 ml of methanol to columns and centrifuge for 1 minute at 200 x g.

Note: Empty tubes to avoid tips of extraction columns from contacting eluates.

- Add 2x 1 ml of distilled or deionized water to columns, centrifuge for 1 minute at 200 x g.
- Proceed with sample application without delay.

SAMPLE APPLICATION

- Load 2x 0.6 ml of EDTA plasma sample onto the correspondingly marked column and centrifuge for 1 minute at 200 x g.

Note: If the EDTA plasma sample is highly viscous, repeat centrifugation for 1 minute at 1000 x g.

WASHING

- Add 1 ml of distilled or deionized water to columns, centrifuge for 1 minute at 500 x g.

Note: If the column is clogged, repeat centrifugation for 1 minute at 1000 x g.

ELUTION OF EXTRACT

- Place each extraction column into a clean correspondingly marked polypropylene tube.
- Add 1 ml of methanol to columns and centrifuge for 1 minute at 200 x g.
- **Collect the eluate and keep it for the evaporation step.**
- Wash column with 2x 1 ml of methanol. Discard the eluted wash solution
- Use column for extracting the next sample (up to 5 times) or store column at 18-28°C and protected from light and dust.

EVAPORATION AND RECONSTITUTION OF EXTRACT

- Evaporate the methanol to dryness using a vacuum concentrator with a cold trap. Alternatively, evaporate the methanol to dryness with a stream of particle free nitrogen.
- Reconstitute the samples with 1.2 ml of Tris buffer and vortex thoroughly.
- Equilibrate the extracts for at least 30 minutes at 2-8°C and vortex again.
- Store reconstituted extracts capped and frozen at -20°C if not assayed immediately.

ASSAY PROCEDURE

Note: Keep assay tubes and reagent vials in an ice bath (0-4°C) during all pipetting steps.

1. Label 8 tubes in duplicate: A to E for the calibrators, NSB for the blank, MB for the maximum binding and T for the total activity tubes. Label additional tubes in duplicate for patient samples and controls.
2. Pipet 600 µl of Tris Buffer into the NSB tubes, and 500 µl of Tris Buffer into the MB tubes.
Pipet 500 µl of the reconstituted Calibrators A to E into the corresponding tubes.
Pipet 500 µl of the extracted patient samples and reconstituted controls into each of the correspondingly marked tubes.
3. Add 100 µl of the Ang II Antiserum to all tubes except the NSB and T tubes. Vortex.
4. Incubate all tubes for 16 hours (+ 4 hours) at 2-8°C.
5. Add 100 µl of ¹²⁵I-Ang II tracer to all tubes. Vortex. Remove the T tubes, they will require no further processing until counting in step 11.
6. Incubate for 6 hours (± 30 minutes) at 2-8°C.
7. Invert the bottle containing the solid phase second antibody several times, add a stir bar and place the bottle on a magnetic stirrer.
While stirring the second antibody suspension continuously, add 100 µl of the suspension to all assay tubes (except the T tubes). Vortex.
8. Incubate for 30 minutes (± 5 minutes) at 2-8°C.
9. Add 1 ml of distilled or deionized water to all tubes (except T tubes).
10. Centrifuge for 5 minutes at 1000 x g and 2-8°C
(*In order to obtain a compact pellet, a centrifugation at 2000-3000 x g during 15 minutes is recommended*).
Aspirate the supernatants (except the T tubes) and retain the precipitates for counting.
11. Count the tubes for 1 minute in a gamma-counter.

RESULTS & Standardization

Record the cpm for all tubes (T, NSB, MB, Calibrators A-E, samples and Controls) and calculate the mean cpm for each pair of tubes. Subtract the mean assay blank (NSB tubes) from the respective mean of each pair of tubes:

$$\text{Net cpm} = \text{cpm}_{\text{Average}} - \text{cpm}_{\text{Average NSB}}$$

Calculate the binding of each pair of tubes as a percent of maximum binding (MB tubes), with the NSB-corrected cpm of the MB tubes taken as 100%.

$$B/B_0(\%) = \text{percentbound} = \frac{\text{netcpm}}{\text{netMBcpm}} \times 100$$

Prepare a lin/log graph paper and plot the percent bound on the vertical axis against the Ang II concentration (pg/ml) on the horizontal axis for each of the Calibrators. Draw the best fitting curve or calculate the standard curve using a four parameter algorithm. Determine Ang II concentrations for the patient samples and Controls from this standard curve. Alternative data reduction methods are equally acceptable.

See Table 11 and Figure 1 for examples of results and standard curve. *These results and standard curve is for*

demonstration purposes only. A standard curve must be generated for each set of samples to be assayed

The Calibrators are calibrated against WHO reference preparation MRC70/302.

QUALITY CONTROL

A thorough understanding of this instruction for use is necessary for the successful use of the product. Reliable results will be obtained only by using precise laboratory techniques (current GLP guidelines) and accurately following this instruction for use.

Ang II is presently not included in any external QC schemes. Therefore, each laboratory should establish its own reference values in addition to the use of an internal control.

The reproducibility of standard curve parameters and control values should be within established limits of laboratory acceptability. The confidence limits for the Controls are lot-specific and printed on the QC Data Sheet added to the kit.

If the precision of the assay does not correlate with the established limits and repetition excludes errors in technique, check the following issues: i) pipetting, temperature controlling and timing devices ii) expiration dates of reagents iii) storage and incubation conditions iv) purity of water.

LIMITATIONS

- Patient samples that are not properly handled, may cause inaccurate Angiotensin II results. EDTA venipuncture tubes are essential to inhibit ACE activity and the subsequent formation of Ang II from Ang I. An angiotensinase inhibitor, such as Bestatin (code: B-A22-BST), is highly recommended to inhibit the transformation of Ang II to Ang III.
- Freezing plasma samples immediately or extracting plasma samples without delay, will preserve the integrity of the Ang II concentration at the time of sampling.
- The interpretation of the results should be examined regarding eventual differences in immunoreactive peptides and biologically active Angiotensin II (2,11).
- Test results should be interpreted in conjunction with information available from clinical assessment of the patient and other diagnostic procedures.
- Some plasma samples may clog extraction columns. Therefore, to clear plasma samples of substances that could interfere with flow rate of the extraction column. It is recommended to centrifuge thawed plasma samples for 1 minute at 10,000 x g and 2-8°C or, alternatively, for 5 minutes at 1000 x g, then proceed to pipet samples onto columns.

PERFORMANCE CHARACTERISTICS

The assay performance characteristics have been validated in duplicates.

Intra-Assay Precision (Within-Run): 8.3%. The intra-assay precision was calculated from the results of 20 pairs of values from four extracted plasma samples in a single run. The values are presented in Table 12.

Inter-Assay Precision (Run-to-Run): 11.5%. Three plasma samples were freshly extracted each time an assay was performed. The inter-assay precision of the samples was then determined from the results of 20 pairs of values in 20 consecutive runs. The values are presented in Table 13.

Detection limit (LoB): 1.0 pg/ml. Twenty Zero Calibrator duplicates (Tris Buffer) were assayed in a single run. The minimal detectable dose of Ang II was calculated by subtracting two standard deviations (SD) from the mean cpm of the maximum binding and intersecting this value with the standard curve obtained in the same run.

Dilution Linearity/Parallelism: 93.5%. Two human plasma samples containing high concentrations of Ang II were diluted

with a pool of Ang II-free human plasma, extracted and subsequently analyzed according to the assay procedure. The values are presented in Table 14.

Spiking Recovery: 99.2%. In the first set of experiments, six plasma samples (4-9) were spiked with different amounts of Ang II, then extracted and analyzed according to the assay procedure. In the second set of experiments, two plasma samples (A and B) were extracted, spiked with 2.5 to 100 pg/ml of Ang II and then assayed according to the assay procedure. The results are presented in Table 15.

Specificity: The cross-reactions of the Ang II antiserum presented in Figure 2 were determined at 50 % binding.

Method Comparison: 46 patients samples were analyzed using the BÜHLMANN Angiotensin II RIA and a commercially available kit for determining plasma renin activity (PRA). A linear regression analysis of the results is presented in Figure 3. The range of values was 0.3 - 81.8 pg/ml for Ang II and 0.1-16.3 ng/ml/hour for PRA, respectively (12).

REFERENCE INTERVALS

In two studies using EDTA plasma samples from 123 apparently healthy donors, the following reference ranges were established (5,12):

These reference ranges should be used as guidelines only. It is recommended that each laboratory establishes its own expected range for its patient population.

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Figure 1: **Example of a Standard Curve**


Table 12: **Intra-Assay Precision (Within-Run)**

| Sample | Mean Value [pg/ml] | SD [pg/ml] | CV [%] |
|----------|--------------------|------------|--------|
| Normal 1 | 5.9 | 0.67 | 11.3 |
| Normal 2 | 5.8 | 0.59 | 10.3 |
| Elevated | 97.6 | 4.24 | 4.3 |
| High | 242.2 | 17.1 | 7.1 |
| Mean | | | 8.3 |

Table 13: **Inter-Assay Precision (Run-to-Run)**

| Sample | Mean Value [pg/ml] | SD [pg/ml] | CV [%] |
|--------------|--------------------|------------|--------|
| Normal | 4.1 | 0.74 | 18.1 |
| Elevated 20 | 24.2 | 2.59 | 10.7 |
| Elevated 150 | 152.9 | 8.84 | 5.8 |
| Mean | | | 11.5 |

Table 14: **Dilution Linearity/Parallelism**

| Sample | Dilution | Observed [pg/ml] | Expected [pg/ml] | O/E [%] |
|--------|----------|------------------|------------------|---------|
| C | 1:1 | 200.9 | --- | --- |
| | 1:2 | 102.0 | 100.5 | 102 |
| | 1:4 | 48.9 | 50.2 | 97 |
| | 1:8 | 21.5 | 25.1 | 86 |
| | 1:16 | 12.0 | 12.6 | 96 |
| | 1:32 | 6.4 | 6.3 | 102 |
| | 1:64 | 3.7 | 3.1 | 117 |
| D | 1:1 | 217.0 | --- | --- |
| | 1:2 | 93.7 | 108.5 | 86 |
| | 1:4 | 49.3 | 54.3 | 91 |
| | 1:8 | 23.2 | 27.1 | 86 |
| | 1:16 | 10.8 | 13.6 | 80 |
| | 1:32 | 6.1 | 6.8 | 90 |
| | 1:64 | 3.0 | 3.4 | 89 |
| Mean | | | | 93.5 |

Table 15: **Spiking Recovery**

| Sample | endogen. Ang II (pg/ml) | Ang II added (pg/ml) | expected value (E) (pg/ml) | observed value (O) (pg/ml) | recovery O/E (%) |
|--------|-------------------------|----------------------|----------------------------|----------------------------|------------------|
| 4 | 5.8 | 10 | 15.8 | 13.4 | 85 |
| 5 | 6.5 | 20 | 26.5 | 24.7 | 94 |
| 6 | 5.4 | 50 | 55.4 | 56.5 | 102 |
| 7 | 6.4 | 85 | 91.4 | 91.2 | 100 |
| 8 | 6.0 | 100 | 106.0 | 105.3 | 99 |
| 9 | 6.4 | 250 | 256.4 | 256.4 | 100 |
| A | 7.1 | 2.5 | 9.6 | 8.9 | 93 |
| | | 5 | 12.1 | 12.0 | 99 |
| | | 10 | 17.1 | 17.8 | 104 |
| | | 25 | 32.1 | 32.3 | 100 |
| B | 12.2 | 50 | 57.1 | 58.0 | 101 |
| | | 5 | 17.2 | 16.3 | 95 |
| | | 10 | 22.2 | 24.0 | 108 |
| | | 20 | 32.2 | 35.1 | 109 |
| | | 50 | 62.2 | 62.1 | 100 |
| | | 100 | 112.2 | 110.4 | 98 |
| Mean | | | | | 99.2 |

Figure 2: **Cross-Reactivity of Ang II Antiserum**


Figure 3: **Method Comparison**


Table 16: **Expected Values**

| | Switzerland (5) | USA (12) |
|-----------------------|-----------------|----------|
| n | 81 | 42 |
| Median Value (pg/ml) | 5.9 | 5.9 |
| S.D. (pg/ml) | 3.41 | 2.82 |
| Lowest Value (pg/ml) | 0.8 | 0.9 |
| Highest Value (pg/ml) | 16.0 | 13.8 |

Table description: cf. "Results" (page 4), "Performance Characteristics" and "Reference Intervals" (page 5).

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APPENDIX II

REFERENCES/ LITERATURREFERENZEN/ REFERENCES/ RIFERIMENTI/ REFERENCIAS

- Emanuel, R.L. *et al.*: Double antibody radioimmunoassay of renin activity and Angiotensin II in human peripheral plasma. *J. Lab. Clin. Med.* **81**, 632-640 (1973).
- Nussberger, J. *et al.*: True versus immunoreactive Angiotensin II in human plasma. *Hypertension* **7** (suppl. 1), 1-7 (1985).
- Mento, P.R. and Wilkes, B.M.: Plasma angiotensins and blood pressure during converting enzyme inhibition. *Hypertension* **6** (suppl. III) 42III-48III (1987).
- Voelker J.R. *et al.*: Improved HPLC-Radioimmunoassay for quantifying Angiotensin II in plasma. *Clin. Chem.* **40**, 1537-1543 (1994).
- Weber, J. and Haldimann, D.: Specific measurement of low Angiotensin II concentrations in plasma. *FASEB J.* **7**, A546 (Abstract 3165) (1993). Poster presented at Exp. Biol. Congress, New Orleans, March 29 - April 1 (1993).
- Unger, T., Gohlke, P.: Tissue renin-angiotensin systems in the heart and vasculature: possible involvement in the cardiovascular actions of converting enzyme inhibitors. *Am. J. Cardiol.* **65**, 31-101 (1990).
- Hidaka H. *et al.*: An improved method for measuring Angiotensin I converting enzyme using a highly sensitive Angiotensin II radioimmunoassay. *Endocrinol. Japon.* **32**, 803-809 (1985).
- Ishida H. *et al.*: Interrelation between the renin-angiotensin system and kallikrein-kinin system in patients with essential hypertension. *Adv. Exp. Med. Biol.* **198**, 329-336 (1986).
- Boer, P. *et al.*: Response of urinary angiotensin to challenges of the renin-angiotensin systems. *Clin. Chim. Acta* **199**, 195-204 (1991).
- Nussberger, J. *et al.*: Octapeptide-specific and sensitive assay for Angiotensin II in plasma. *Int. J. Environ. Anal. Chem.* **25**, 257-268 (1986).
- Semple, P.F. *et al.*: Angiotensin II and its heptapeptide (2-8), hexapeptide (3-8) and pentapeptide (4-8) metabolites in arterial and venous blood of man. *Circ. Res.* **39**, 671-678 (1976).
- Sica D. *et al.*: A novel analytical method for the determination of Angiotensin II. *Clin. Chem.* **42**, S231 (Abstract 565) (1996). Poster presented at the 48th Annual Meeting of the Am Assoc Clin Chem, Chicago, July 28 - August 1 (1996).

| PIPETTING PROTOCOL ANGIOTENSIN II | | | | | | | | |
|-----------------------------------|------------------|--------------------------------|----------------|--|-------------|--|----------------------|--|
| Polystyrene Tubes in Duplicate | Tris Buffer (μl) | Standard, Control, Sample (μl) | Antiserum (μl) | | Tracer (μl) | | Second Antibody (μl) | |
| Total | | | | | 100 | | | vortex and incubate at 2-8°C for 30 min (± 5 min) |
| NSB | 600 | | | | 100 | | 100 | |
| MB | 500 | | 100 | | 100 | | 100 | add 1 ml of deionized water and centrifuge for 5 min at 2-8°C and 1000 x g |
| Std 2 pg/ml | | 500 | 100 | | 100 | | 100 | |
| Std 5 pg/ml | | 500 | 100 | | 100 | | 100 | vortex and incubate at 2-8°C for 6 hours (± 30 min) |
| Std 20 pg/ml | | 500 | 100 | | 100 | | 100 | |
| Std 100 pg/ml | | 500 | 100 | | 100 | | 100 | |
| Std 500 pg/ml | | 500 | 100 | | 100 | | 100 | aspirate supernatant (except T tubes) and count for 1 min |
| Normal Control | | 500 | 100 | | 100 | | 100 | |
| Elevated Control | | 500 | 100 | | 100 | | 100 | |
| Sample | | 500 | 100 | | 100 | | 100 | |

Table 17

APPENDIX IV
SYMBOLS/ SYMBOLE/ SYMBOLES/ SIMBOLI/ SIMBOLOS

| Symbol | Explanation |
|---|--|
|  | Use By Verwendbar bis Utiliser jusqu'au Utilizzare entro Fecha de caducidad |
| REF | Catalogue number Bestellnummer Référence du catalogue Numero di catalogo Número de catálogo |
| LOT | Batch code Chargenbezeichnung Code du lot Codice del lotto Codigo de lote |
|  | Temperature limitation Zulässiger Temperaturbereich Limites de température Limiti di temperatura Limite de temperatura |
|  | Consult Instructions for Use- Gebrauchsanweisung beachten Consulter le mode d'emploi Consultare le istruzioni per l'uso Consulte las instrucciones de uso |
|  | Contains sufficient for <n> tests Ausreichend für „n“ Ansätze Contenu suffisant pour „n“ tests Contenuto sufficiente per „n“ saggi Contenido suficiente para <n> ensayos |
|  | Radioactive Material Radioaktives Material Matériel radioactif Materiale radioattivo Material radiactivo |

| Symbol | Explanation |
|----------------------|--|
| AEC | Extraction Columns Extraktionssäulen Colonne d'extraction Colonne d' estrazione Columnas de extracción |
| BUF TRIS | Tris Buffer Tris-Puffer Tampon Tris Tampone Tris Tampón Tris |
| Ab | Antiserum Antiserum Antisérum Antisiero Antisuero |
| TR | Tracer Tracer Traceur Elemento tracciante Trazador |
| CAL A - CAL E | Calibrator A - E Kalibrator A - E Calibreur A - E Calibratore A - E Calibrador A - E |
| CONTROL N | Control Normal Normalkontrolle Contrôle normal Controllo normale Control normal |
| CONTROL H | Control High Kontrolle hoch Contrôle élevé Controllo alto Control alto |
| Ab2 | 2 nd Antibody 2. Antikörper 2 ^{ème} Anticorps Secondo anticorpo Segundo anticuerpo |

Printing Date
2013-01-14