Cystatin C

PETIA

Particle enhanced Turbidimetric Immunoassay

KK-CYC
KK-CYC-6*

*Kit contains pre-diluted Calibrators for the use in the clinical chemistry analysers of Olympus (AU) and Ortho (Vitros 5.1)

Revision date: 2008-11-13

ENGLISH

INTENDED USE

The Cystatin C Immunoassay is an *in-vitro* diagnostic test for quantitative determination of cystatin C in human serum and plasma. The assay can be applied on any clinical analyzer according to the specific application note. The measurement of cystatin C is used in the diagnosis and treatment of renal diseases.

SUMMARY AND EXPLANATION OF TEST

The non-glycosylated basic protein, cystatin C (molecular weight 13.2 kDa), is produced at a constant rate in nearly every nucleated cell in the human body [1]. It is freely filtered through a normal Glomerular membrane, and is then reabsorbed and almost entirely catabolized in the proximal tubules. Hence, the cystatin C concentration in human blood is closely related to glomerular filtration rate (GFR) [2]. A reduction in the GFR causes a rise in the concentration of cystatin C. Cystatin C concentration has not been shown to be significantly influenced by other factors such as muscular mass, inflammatory diseases, sex, age or diet [2, 3, 4].

ASSAY PRINCIPLE

Serum or plasma sample from human is mixed with cystatin C immunoparticles. Cystatin C from the sample and anti cystatin C from the immunoparticles aggregate. The complex particles created absorb light, and by turbidimetry the absorption is related to cystatin C concentration via interpolation on an established standard calibration curve.

REAGENTS SUPPLIED AND PREPARATION

Reagents	Quantities	Code	Reconstitution
Reaction Buffer (R1) a MOPS [3-(N-Morpholino)-propanesulfonic acid] buffered saline, preserved with sodium azides (0.09 % (w/v))	1 vial 55 ml	B-CYC-R1	Ready to use
Immunoparticles (R2) purified immunoglobulin fraction directed against Cystatin C. Covalently attached to uniform polystyrene particles. preserved with 0.09 % (w/v) sodium azide and antibiotics	1 vial 10 ml	B-CYC-R2	Ready to use
Calibrator¹⁾ delipidated human serum pool spiked with human cystatin C. Antibiotics are used as preservative.	1 vial 1 ml	B-CYC-CA	Ready to use
Calibrator pre-diluted²⁾ Pre-diluted Calibrators A-F	6 x 1 vial 1 ml	B-CYC-CASET	Ready to use
Controls low/high¹⁾ delipidated human serum pool spiked with human cystatin C. Antibiotics are used as preservative.	2 vials 1 ml	B-CYC-CONSET	Ready to use

Table 1

¹⁾ Content of KK-CYC.

Calibrator values are lot dependent and a new calibration must be performed whenever a new calibration lot is used. The calibrator's assigned values are given on the analytical value sheet provided with the kit.

²⁾ Content of KK-CYC-6

The pre-diluted Calibrator set (B-CYC-CASET) will be needed for the application on the Olympus AU and the Vitros 5,1 Analyser from Ortho. Applying the Cystatin C assay on these analyzers please order the kit KK-CYC-6

The reagents last for at least 220 tests dependent on the analyzer specific instrument settings (see section Materials required)

STORAGE AND SHELF LIFE OF REAGENTS

Unopened Reagents	
Store at 2-8°C until expiration date. Do not use past expiration date.	
Opened Reagents	
Reaction Buffer	Stable for 9 weeks at 2-8°C On board stability: 4 weeks at 2-8°C
Immunoparticles	Reagents should be stored capped at 2-8°C when not in use.
Calibrator(s)	Store at 2-8°C. stable until expiration date.
Controls	

Table 2

WARNINGS AND PRECAUTIONS

- This test is for *in vitro* use only, and must be handled by qualified personnel.
- Reagents contain antibiotics and must be handled with due caution.
- Reagents contain sodium azide preservative and must be handled with due caution: Do not ingest or allow contact to skin or mucous membranes. The sodium azide concentration of this product is not characterised as dangerous. Although, accumulated NaN₃ in lead and copper pipes may cause generation of explosive metal azides. To prevent this, rinse thoroughly if discarded into the drain.
- The immunoparticles contain substances of animal origin. Disposal of any discarded materials should be in accordance with local requirements.
- Serum used in the manufacture of calibrators and controls was tested for hepatitis HBsAG, anti-HCV, anti-HIV1 and anti-HIV2 and found to be negative. Nevertheless, the materials contain substances of human and animal origin and must be handled with due care. Disposal of any discarded materials should be in accordance with local requirements.

MATERIALS REQUIRED BUT NOT PROVIDED

- Analyzer specific reagent cartridge for R1 and R2 as well as for Calibrator and Controls
- instrument specific settings (parameter lists) are available for the following analyser upon request:

Hitachi 911	Advia 1650/1800
Hitachi 917	Olympus AU*
Roche Modular P	Ortho Vitros 5,1*
Roche Cobas c501	Abbott Architect

A General Application Note is also available

*The instrument specific settings for Olympus AU and Ortho Vitros 5,1 require pre-diluted calibrators which will be provided by BÜHLMANN AG with the Cystatin C kit Order code: KK-CYC-6.

SPECIMEN COLLECTION AND HANDLING

Required sample material is human serum or EDTA/Heparinized plasma. It is recommended to analyse fresh samples. Mix samples well before analysing. Sample stability testing shows that cystatin C in serum and plasma samples are stable for 14 days at room temperature (8-25°C) and for 21 days if stored at 2-8°C, otherwise store frozen at ≤-20°C [8].

ASSAY PROCEDURE

Application Notes/Assay Installation

An Application Note is available for the analyzer mentioned above and will be provided upon request. Instrument set-up, maintenance, operation and precautions must be handled in accordance with the instrument manuals.

Reagent Preparation

Reagents are supplied ready to use. Transfer the content of Reagent 1 (Reaction buffer) and Reagent 2 (Immunoparticles) into appropriate compartments of the User Defined Cartridge as described in the Application Note. Use care to avoid contamination. Reagents should be stored capped at 2-8°C when not in use.

Establishment of the Calibration Curve

Use BÜHLMANN Cystatin C Calibrator (B-CYC-CA or B-CYC-CASET) to establish a 6-point standard curve as defined in the instrument manuals. Calibrator values are lot dependent and a new calibration must be performed whenever a new calibration lot is used. The calibrator's assigned values are given on the analytical value sheet provided with the kit.

QC Controls

The controls (low and high) must be assayed each day before any samples are assayed in order to validate the calibration curve. The controls have an assigned value range that must be met before measuring samples. The assigned values are given in the lot specific data sheet included with the BÜHLMANN Cystatin C kit. If the control values are not valid, repeat the control measurements. If the calibration cannot be performed without error, or valid control values cannot be reproduced, contact the instrument manufacturer for support.

Measuring Patient Samples

When a valid calibration has been performed and the control values are within the valid range, serum or plasma samples may be measured. Check that minimum volume of sample is present and assay the samples according to the instructions given in the Instrument Manuals.

RESULTS

The results are presented in mg/L.

The GFR values can be calculated using the **GFR Calculator v2.0** on the website: www.gfrcalc.com

Measuring range

The measuring range of Cystatin C for the assay is approximately 0.4 – 8.0 mg/L. The exact range is dependent on the calibrator value, which is lot specific.

LIMITATIONS

The materials shall not be used past expiration date. Interference studies are a part of instrument validation. For a variety of analysers, no interference is detected for this product with triglycerides, haemoglobin, Intralipid or bilirubin at tested concentrations. There is no interference detected with the drugs tested on recommendation from Sonntag and Scholer [5]. There is no RF interference present because the antibodies are made from an avian source (chicken) [6].

REFERENCE INTERVALS

The CLSI Guideline, C28-A2; *How to Define and Determine Reference Intervals in the Clinical Laboratory; Approved Guideline Second edition* was followed to determine the transferability of the reference interval. For Architect® ci8200 the reference interval is 0.53-1.01mg/L. For any other instrument the reference interval should be in

accordance with the interval obtained on the Architect® ci8200. It is recommended that every laboratory determine a local reference interval since values may vary depending on the population tested.

PERFORMANCE CHARACTERISTICS

All results refer to the validation of BÜHLMANN Cystatin C Immunoassay on Abbott Architect® ci8200.

Detection Limit

Limit of Quantification is defined as: *The lowest actual amount of an analyte that can be reliably detected and at which the total error meets the requirements for accuracy.*

Limit of Detection is defined as: *The smallest amount of an analyte that the method can reliably detect to determine presence or absence of the analyte.* BÜHLMANN Cystatin C Immunoassay has a Limit of Detection of 0.031mg/L and a Limit of Quantification of 0.33 mg/L.

Precision

The BÜHLMANN Cystatin C Immunoassay was used in a 5 day precision study designed in accordance with CLSI protocol EP5-A. 4 serum pools and 2 control levels were measured on the Architect® ci8200.

Sample	Mean value (mg/L)	Within- Run CV (%)	Total CV (%)
Low human serum pool	0,69	1,51	2,81
High human serum pool	5,71	1,12	4,18
Medium human serum pool	3,38	1,93	3,75
Medium human serum pool	1,35	0,97	2,67
Low Control	0,88	1,42	3,72
High Control	3,58	0,71	1,41

Linearity

BÜHLMANN Cystatin C Immunoassay on Architect® ci8200 is linear in the range of 0.3 – 8.8 mg/L. Concentrations outside this range have not been tested.

Analytical Recovery

For the BÜHLMANN Cystatin C Immunoassay on Architect® ci8200 a recovery of 100-107% can be expected.

Interference

No interference is detected with Triglycerides (12.5 mmol/mL) Haemoglobin (8.0 g/L), Intralipid (11.0 g/L) or Bilirubin (420 mg/L) on Architect® ci8200 .

There is no interference detected with the drugs tested on recommendation from Sonntag and Scholer [5]. The interference study was designed in accordance with the protocol EP7-A from CLSI [7]. There is no RF interference present in the BÜHLMANN Cystatin C Immunoassay because the antibodies are made from avian (chicken) [6].

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