



ACE kinetic

Procedure

KK-ACK


A Commitment to Diagnostics

INTENDED USE

BÜHLMANN ACE kinetic is an *in vitro* diagnostic bio-chemical assay for the quantitative determination of angiotensin converting enzyme (ACE) activity in serum samples. The assay aids the assessment of disease activity in patients with sarcoidosis in conjunction with other clinical and laboratory findings.

PRINCIPLE OF THE ASSAY

The assay is a quantitative enzymatic test which can be easily applied on clinical chemistry analyzers or run by manual method. ACE catalyzes the conversion of angiotensin I to angiotensin II. The enzyme also mediates the cleavage of the synthetic substrate FAPGG (N-[3-(2-furyl)acryloyl]-L-phenylalanyl-L-gly-cyl-L-glycine) into the amino acid derivative FAP and the dipeptide GG. The linear kinetic of this cleavage reaction is measured by recording the decrease in absorbance at 340 nm. The final ACE activity in U/L in the patient sample is determined using a calibration curve generated from the measured calibrator value (Ronca-Testoni, Clin Chem 1983; Bénétteau, Clin Chem 1986).



ASSAY PROCEDURE

Application notes / assay installation

Assay procedures for the ACE kinetic are established on several clinical chemistry analyzers. Validated application notes describing installation and analysis on specific instruments are available from BÜHLMANN upon request. Corresponding instrument manuals must be considered for instrument setup, maintenance, operation and precautions.

Validated Applications

Roche cobas®	c501/502 c701/702 c303 (Pure) c503 (Pro)
Abbott	Alinity c Architect c-series
Beckman	AU480/AU680 AU5800/DxC700AU
Siemens	Atellica CH930 Advia 2400 Dimension Vista Dimension EXL
The Binding Site	Optilite
ThermoFisher Scientific	Indiko
IDS	iSYS
Mindray	BS480

For other clinical chemistry analyzers please contact info@buhlmannlabs.com.

Manual procedure on microtiter plate is possible applying an MP Reader with 37°C incubation and plate shaking option and optical filter at 340nm and 415 nm.

Pre-Analytcs

Sample required:	~200 µL serum
Sample collection	Gel separator tubes (SST) are recommended Optionally, Li-Heparin and Citrate Plasma can be used* EDTA Plasma inhibits ACE activity Serum collection tubes without anti-coagulants
Sample storage:	at 2-8°C up to 30 days at -20°C at least 6 months

*Plasma samples not validated. Plasma testing is at discretion of the laboratory.

Special Equipment

Open clinical chemistry analyzer: validated applications are available

For manual procedure: kinetic spectrophotometer with 340 and 415 nm filter and incubation chamber at 37°C

Kit Components

ACE kinetic is available in different package sizes.

	KK-ACK	KK-ACK4	KK-ACK2	KK-ACKX
Tests	100	400	2 x 50	1200
Substrate	1 x 26 mL	4 x 26 mL	2 x 13 mL	3 x 100 mL
Calibrator	1 x 2 mL	2 x 2 mL	2 x 2 mL	3 x 2 mL
Controls normal/high	1 x 2 mL	2 x 2 mL	2 x 2 mL	3 x 2 mL





ACE kinetic

Characteristics

KK-ACK

A Commitment to Diagnostics

Assay Performance

Data obtained on the Roche cobas® c501. Refer to specific application notes for performance on other clinical chemistry analyzers.

Reproducibility **6.3-9.1% CV**
3 instruments/lots x 5 days x 5 replicates (EP05-A3)
On cobas c501, c701 and AU480

Repeatability **0.8-3.0% CV**
Within-Laboratory Precision 1.7-3.7% CV
20 days x 2 runs x 2 replicates (EP05-A3)

Accuracy/Recovery **92.0 – 112.8%**
6 samples spiked with 20.5 U/L (10% volume) run in 4 replicates

Sample carry-over
No statistical significant carry-over (EP10-A2)

Limit of Blank (LoB) **4.3 U/L**

Limit of Detection (LoD) **6.8 U/L**

Limit of Quantification (LoQ) **11.3 U/L**
Classical approach EP17-A2; LoQ n=60; <20% CV

Linearity range **4.3-535 U/L**
Samples >150 U/L automatically re-run with reduced volume; acceptance ±4 U/L or ±10% (EP06-A)

Security zone **up to 541 U/L**
No limiting effects observed

INTERFERING SUBSTANCES

Susceptibility to interfering substances assessed according EP07-A2. Bias exceeding 20% considered as interference.

Oral pharmaceuticals

No interference detected:

- Aspirin 0.65 mg/mL
- Azathioprine 3.0 µg/mL
- Chlorambucil 7.2 µg/mL
- Cyclophosphamide 0.375 mg/mL
- Eprosartan 0.36 mg/mL
- Hydroxychloroquine up to 0.06 mg/mL
- Ibuprofen 0.5 mg/mL
- Losartan 0.09 mg/mL
- Methotrexate 2.0 µg/mL
- Prednisone 0.3 µg/mL

Serum indices

Interference detected at concentrations:

- triglycerides 2.24 mg/mL
- conjugated bilirubin 0.06 mg/mL
- unconjugated bilirubin 0.047 mg/mL
- hemoglobin 1.19 mg/mL

No triglycerides interference observed after short centrifugation (10 min / 12'000 x g) and separation of lipid-containing supernatant.

This document is for information purpose only, before performing the assay please carefully refer/read the respective IFU available (<https://www.buhlmannlabs.com/products-solutions/clinical-chemistry/ace-kinetic>).

REFERENCE INTERVALS

Adults

2.5th – 97.5th percentile from healthy participants in three independent studies in Switzerland (n=80, age: 20 – 70), Germany (n=159, age: 18 – 64, ref. 3) and USA (n=327, age 16 – 77):

20 – 70 U/L

Children

2.5th – 97.5th percentile from healthy pediatric participants in a single study in Germany (n=84, age: 0.5 – 18):

33 – 112 U/L

Plasma samples*

Samples from healthy blood donors collected into lithium-heparin and citrate tubes compared to serum samples from the same donors:

Li-Hep Plasma (n=38) $y=0.9x + 2.5$; $r=0.975$

Mean bias: -1.1%

Citrate Plasma (n=44) $y=0.8x + 1.7$; $r=0.990$

Mean bias: -10.8%

*Plasma samples not validated. Plasma testing is at discretion of the laboratory.

Ordering codes:

- KK-ACK 100 tests
- KK-ACK2 2x50 tests
- KK-ACK4 400 tests
- KK-ACKX 1200 tests

ACE Kinetic is FDA 510(k) Exempt. For *in vitro* Diagnostic Use.

CE 0123



BÜHLMANN Laboratories AG
Switzerland
info@buhlmannlabs.ch
www.buhlmannlabs.com
Phone +41 61 487 12 12

BUHLMANN Diagnostics Corp
USA
info@buhlmannlabs.com
www.buhlmannlabs.com
Phone 844 300 9799