# CIC-C1Q

## Circulating Immune Complexes

### ELISA

EK-CIC      96 wells

Revision date: 2012-07-23

**INTENDED USE**

The BÜHLMANN CIC-C1Q EIA is intended for the quantitative *in vitro* diagnostic determination of circulating immune complexes (CIC) in serum and plasma of patients with various autoimmune and other CIC-related diseases (1-4).

**PRINCIPLE OF THE ASSAY**

Circulating immune complexes from patient sera, plasma or Calibrators and controls are incubated with human C1q adsorbed onto microtiter wells. After a washing step, an alkaline phosphatase (aPase) labeled conjugate (Protein A) is added, which binds to the Fc region of human IgG. After another washing step, the enzyme substrate (*para*-nitrophenyl-phosphate or pNPP) is pipetted followed by a stopping step (5-7). The absorption is measured at 405 nm.

**REAGENTS SUPPLIED AND PREPARATION**

Reagents	Quantity	Code	Reconstitution
<b>Microtiter Plate</b> Precoated with human C1q	12 x 8-well	B-CIC-MP	Ready to use
<b>Plate Sealer</b>	3 pieces		
<b>Wash Buffer Concentrate (10x)</b> with preservatives	1 bottle 100 ml	B-CIC-WB	Dilute with 900 ml of deionized water
<b>Incubation Buffer</b> with preservatives	1 bottle 100 ml	B-CIC-IB	Ready to use
<b>Calibrators A to E<sup>1)</sup></b> Aggregated-IgG in a buffer matrix with preservatives	5 vials 1 ml	B-CIC-CASET	Ready to use
<b>Control Low / High<sup>2)</sup></b> Aggregated-IgG in a buffer matrix with preservatives	2 vials 1 ml	B-CIC-CONSET	Ready to use
<b>Enzyme Label</b> Proteine A conjugated aPase in a protein-based buffer matrix with preservatives	1 vial 11 ml	B-CIC-EL	Ready to use
<b>pNPP Substrate</b> with stabilizing pellets	1 vial 11 ml	B-PNPP	Ready to use
<b>Stop Solution</b> 1 N sodium hydroxide	1 vial 11 ml	B-NAOH	Ready to use <b>Corrosive agent</b>

Table 1

<sup>1)</sup> The Calibrators A ,B, C, D and E contain 50, 20, 10, 4 and 1 µg Equivalents of aggregated-IgG per ml (µg Eq/ml), respectively.

<sup>2)</sup> The controls contain lot-specific amounts of Equivalents of aggregated-IgG per ml (µg Eq/ml). Refer to the QC data sheet added to the kit for exact concentrations.

**STORAGE AND SHELF LIFE OF REAGENTS**

Unopened Reagents	
Store at 2-8°C <b>except Calibrators and Controls</b> which must be stored at <b>-20°C</b> or below. Do not use past kit expiration date.	
Opened / Reconstituted Reagents	
Microtiter Plate	Return unused strips immediately to the foil pouch containing the desiccant packs and reseal along the entire edge of zip-seal. Store until expiration date at 2-8°C.
Wash Buffer	Store at 2-8°C until expiration date.
Calibrators	Stable at -20°C until expiration date. For storage shorter than 3 months Calibrators and Controls can be stored at 2-8°C
Controls	
Incubation Buffer	Store at 2-8°C until expiration date.
Enzyme Label	
Substrate Solution	
Stop Solution	

Table 2

**SAFETY PRECAUTIONS**

- The microtiter-plate, calibrators and controls of this test contain components of human origin. Material of human may contain infectious agents e.g. Hepatitis B virus and HIV.
- **SUBSTRATE AND STOP SOLUTION:** Avoid contact with Eyes, skin and cloths Wear suitable protective clothing, gloves and eye protection. After contact with eyes or skin, wash immediately with plenty of water
- Unused solutions should be disposed of according to local State and Federal regulations.

**TECHNICAL PRECAUTIONS**

- Components must not be used after the expiry date printed on the labels. Do not mix different lots of reagents.
- Avoid contamination of reagents.
- Microwells cannot be re-used.
- It is important reading through the Instructions for Use prior to commencing the test. Reliable results will be obtained only if this Instruction for Use is followed accurately.
- The substrate solution is ready to use. **Do not vortex or try to homogenize the stabilizing pellets prior to use.**

**MATERIALS REQUIRED BUT NOT PROVIDED**

- Precision pipettes with disposable tips for 20 µl, 100 µl and 1000 µl.
- Disposable polystyrene or polypropylene tubes for the preparation of sample dilutions.
- 1000 ml cylinder for the dilution of the Wash Buffer concentrate.
- Microtiter plate washer or squeeze bottle for Wash Buffer.
- Blotting paper.
- Microtiter plate rotator.
- Microtiter plate reader for measurement of absorbance at 405 nm.

**SPECIMEN COLLECTION AND STORAGE**

The procedure calls for 50 µl of blood per duplicate determination. Collect blood into plain tubes, avoid hemolysis, mix by inverting sample tube several times and leave to clot for 45 minutes at room temperature (18-28°C) protected from light. Centrifuge at 1800 x g for 15 minutes at room temperature (18-28°C) and collect the serum. Samples may be stored at 2-8 °C for up to 7 days. If samples are to be stored for a longer period of time, they should be kept at ≤-20°C.

**ASSAY PROCEDURE**

**Important: Allow all reagents to 18-28°C prior to use.**

1. Dilute all patient samples 1:50 with Incubation Buffer (e.g. 20 µl of serum or plasma + 980 µl of Incubation Buffer) and mix well. Allow diluted samples to stand for 15 minutes at 18-28°C prior to pipetting in step 4c.
2. Prepare a plate with sufficient strips to test the required number of calibrators, controls and samples. Remove excess strips from the holder and re-seal them in the foil pouch together with the desiccant packs **without delay**. Store refrigerated.
3. Wash the coated wells twice using at least 300 µl of Wash Buffer per well. Empty the wells and strike the plate firmly onto blotting paper.
- 4a. Pipet 100 µl of Incubation Buffer in duplicate into wells A1+A2

- Pipet 100 µl of Calibrator A in duplicate into wells B1+B2  
 Pipet 100 µl of Calibrator B in duplicate into wells C1+C2.  
 Pipet 100 µl of Calibrator C in duplicate into wells D1+D2.  
 Pipet 100 µl of Calibrator D in duplicate into wells E1+E2.  
 Pipet 100 µl of Calibrator E in duplicate into wells F1+F2.
- 4b. Pipet 100 µl of the Low Control in duplicate into wells G1+G2.  
 Pipet 100 µl of the High Control in duplicate into wells H1+H2.
  - 4c. Pipet 100 µl of each diluted sample (1:50) in duplicate into the subsequent wells.
  5. Cover the plate with a plate sealer, place the plate on a plate rotator set at 400-600 rpm and incubate for 1 hour  $\pm$  5 minutes at 18-28°C.
  6. Remove and discard the plate sealer. Empty the wells and wash three times using at least 300 µl of Wash Buffer per well. Empty the wells and strike the plate firmly onto blotting paper.
  7. Pipet 100 µl of Enzyme Label to all wells.
  8. Cover the plate with a plate sealer, place the plate on a plate rotator set at 400-600 rpm and incubate for 30 minutes  $\pm$  5 minutes at 18-28°C.
  9. Remove and discard the plate sealer. Empty the wells and wash three times using at least 300 µl of Wash Buffer per well. Empty the wells and strike the plate firmly onto blotting paper.
- Important: Allow the pNPP substrate solution to reach 18-28°C.**
10. Pipet 100 µl of the pNPP Substrate Solution to all wells.
  11. Cover the plate with a plate sealer, place the plate on a plate rotator set at 400-600 rpm, protect the plate from direct light and incubate for 30  $\pm$  5 minutes at 18-28 °C.
  12. Pipet 100 µl of Stop Solution to all wells. Remove air bubbles with a pipette tip. Proceed to step 13 within 30 minutes.
  13. Read the absorbance at 405 nm in a microtiter plate reader.

## RESULTS AND CALCULATION

### Standard Curve

- Record the absorbance at 405 nm for each calibrator and blank well.
- Average the duplicate values, subtract the average of the blank wells and record averages (= corrected average absorbance).
- Plot the absorbance (vertical axis) versus the CIC concentration in µg Eq/ml of the Calibrators (horizontal axis) using log/log graph paper.
- Draw the best fitting curve or calculate the standard curve using a four parameter algorithm.

### Samples and Controls

- Record the absorbance at 405 nm for each sample and control well.
- Average the duplicate values, subtract the average of the blank wells and record the averages (=corrected average absorbance).
- Locate the corrected absorbance value of the sample on the vertical axis, draw a horizontal line intersecting the standard curve and read the CIC concentration (µg Eq/ml) from the horizontal axis.

See Table 11 and Figure 1 for an example of results and standard curve. *These results and standard curves are for demonstration purposes only. A standard curve must be generated for each set of samples to be assayed.*

## QUALITY CONTROL

A thorough understanding of this package insert is necessary for the successful use of the product. Reliable

results will be obtained only by using precise laboratory techniques (current GLP guidelines) and accurately following this package insert.

Since there is no control serum for CIC commercially available, we recommend using a positive serum pool for internal quality controls.

All controls must fall within established confidence limits. The confidence limits for the Controls are lot-specific and printed on the additional data sheet.

The reproducibility of standard curve parameters and control values should be within established limits of laboratory acceptability. If the precision of the assay does not correlate with the established limits and repetition excludes errors in technique, check the following issues: i) pipetting, temperature controlling and timing devices ii) ELISA reader settings iii) expiration dates of reagents iv) storage and incubation conditions v) pNPP Substrate Solution should be colorless vi) purity of water.

## PERFORMANCE CHARACTERISTICS

**Intra-Assay Precision (Within-Run): 3.6%.** The intra-assay precision was calculated from the results of 20 pairs of values obtained in a single run. The values are presented in µg Eq/ml (see Table 12)

**Inter-Assay Precision (Run-to-Run): 11.3%.** The inter-assay precision was calculated from the results of 3 pairs of values obtained in 20 different runs. The values are presented in µg Eq/ml (see Table 13).

**Dilution Linearity/Parallelism: 104.8%.** Human serum samples containing low and high titers of CIC were diluted with Incubation Buffer and subsequently assayed according to the assay procedure. The values are presented in µg Eq/ml (see Table 14).

**Plasma/Serum comparison:** Human plasma and serum samples from 47 donors were diluted (1:50) with Incubation Buffer and subsequently assayed according to the assay procedure. All negative (n=27) and all borderline (n=10) serum samples were negative or borderline, respectively when measured in EDTA plasma. The same was true for heparin plasma (n=20). Ten positive serum samples showed a very high correlation (r=0.978) to the plasma samples of the same donor. Therefore CIC can be tested either in serum or plasma samples.

**Standardization:** The Calibrators of the BÜHLMANN CIC-C1Q EIA Kit were calibrated against a reference preparation of aggregated human IgG, which was prepared under the auspices of WHO (Geneva Section).

### Detection Limit:

**Limit of Blank (LoB): 0.58 µg Eq/ml.** Twenty duplicates of Incubation Buffer (Reagent Blank) were assayed in a single run. Mean and standard deviation were calculated for the absorbance values. The minimum detectable dose of CIC complexes was calculated to be 0.58 µg Eq/ml by adding two standard deviations to the mean absorbance of the Reagent Blank (Incubation Buffer) and intersecting this value with the standard curve obtained in the same run.

**Limit of Quantification (LoQ): <1 µg Eq/ml.** Twenty duplicates of low CIC concentrations were assayed in a single run. The limit of quantification was calculated to be <1 µg Eq/ml (cut-off intra-assay CV = 10%).

**Method of Comparison:** A sample comparison (n=64) of this assay and another FDA cleared commercial method showed the following correlation:

Relative Sensitivity	24/25 = 96%
Relative Specificity	38/39 = 97.4%
Relative Accuracy	62/64 = 96.9%

---

## CUT-OFF VALUES (EXPECTED VALUES)

The frequency of CIC in normal human sera or plasma was determined using blood samples from asymptomatic volunteer blood donors (adult men and women at the age of 18 to 70 years). 192 samples were assayed according to the assay procedure and the results shown in Table 15, expressed in  $\mu\text{g}$  Equivalent per ml ( $\mu\text{g}$  Eq/ml), were obtained.

**Proposed Cut-off Titer:** A total of 10 strongly elevated values ( $> \text{mean}+3\text{SD}$ ) from the results of the 192 apparently healthy normal blood donors were eliminated in three iterative cycles. This resulted in a theoretical **cut-off value** ( $\text{mean} + 2 \text{SD}$ ) of **3.2  $\mu\text{g}$  Eq/ml**.

**Borderline and positive:** For practical reasons we recommend that values between **3.2 and 5.0  $\mu\text{g}$  Eq/ml** should be regarded as **borderline results** (grey zone).

Values above **5.0  $\mu\text{g}$  Eq/ml** should be regarded as **positive**.

**General limitations of CIC normal ranges:** elevated ranges of CIC are known in up to 10% of normal blood donors without any clinical manifestations (true negatives outside the normal range) (8,9). This stands in good agreement with our results: 6 (3.1 %) of 192 normal donor samples were positive in the CIC ELISA and 14 (7.3 %) were within the grey zone.

---

## PERFORMANCE LIMITATIONS

- **Interferences:** Lipemic, hemolytic and icteric samples should not be used in this assay. Lipemic samples can be avoided by asking patients to fast for at least 12 hours prior to the sample being taken.
- **Application:** The reagents supplied with this kit are optimized to measure circulating immune complexes in human serum or plasma.
- **Assessment:** The amount of circulating immune complexes should be used as supplementary data available to the physician in establishing a diagnosis.
- **Results.** Serum or plasma titer values depend on the assay method and, in particular, on the specificity and the cut-off values established with the assay method. Titer values obtained with different assay methods cannot be compared directly.
- **Hook effect:** No hook effect was observed up to 50'000  $\mu\text{g}$  Eq/ml.
- **Expected values:** The proposed cut off titers should be used as guidelines only. It is recommended that each laboratory establishes its own reference intervals.

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
Figure 1: **Example of Standard Curve**Table 12: **Intra-Assay Precision (Within-Run)**

Sample #	Mean	SD	CV (%)
1	11.6	0.71	6.13
2	19.0	0.70	3.66
3	19.2	0.38	1.98
4	22.5	0.59	2.65
<b>Mean:</b>			<b>3.61</b>

Table 13: **Inter-Assay Precision (Run-to-Run)**

Sample	Mean	SD	CV (%)
5	9.2	1.1	12.0
6	14.6	1.8	12.6
7	29.4	2.7	9.2
<b>Mean:</b>			<b>11.3 %</b>

Table 14: **Dilution Linearity/Parallelism**

Sample #	Dilution	Observed	Expected	O/E (%)
13	1/ 50 100%	5.27		
	1/100 50%	3.19	2.64	121
	1/200 25%	1.49	1.32	113
	1/400 12.5%	0.84	0.66	127
14	1/ 50 100%	9.86		
	1/100 50%	4.83	4.93	98
	1/200 25%	2.41	2.47	98
	1/400 12.5%	1.13	1.23	91
15	1/ 50 100%	21.84		
	1/100 50%	11.96	10.92	110
	1/200 25%	5.15	5.46	94
	1/400 12.5%	2.49	2.73	91
<b>Mean:</b>				<b>104.8</b>

Table 15: **Cut-Off Values**

	Before Elimination	After Elimination
Total (n)	192	182
Median Value (µg Eq/ml)	1.25	1.19
Mean Value (µg Eq/ml)	1.68	1.50
SD (µg Eq/ml)	1.34	0.83
Lowest Value (µg Eq/ml)	0.52	0.52
Highest Value (µg Eq/ml)	8.9	4.14

**Table description:** cf. "Results and Calculation" (page 3), "Performance Characteristics" (page 3) and "Cut-Off Values" (page 4).

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1. V
2. an Hoeyveld E. & Bossuyt X.: *Evaluation of seven commercial ELISA kits compared with the C1q solid-phase binding RIA for detection of circulating immune complexes.* Clin Chem. **46**, 283-285 (2000).
3. Stanilova S.A. & Slavov E.S.: *Comparative Study of circulating immune complexes quantity detection by three assays – CIF-ELISA, C1q-ELISA and anti-C3 ELISA.* J Immunol Meth **253**, 13-21 (2001).
4. Nydegger, U.E.: *A place for soluble immune complexes in clinical immunology.* Immunol Today **6** (3) 80-82 (1985).
5. Hay, F.C. et al.: *Routine assay for the detection of immune complexes of known immunoglobulin class using solid phase C1q.* Clin. Exp. Immunol. **24**, 396-400 (1976).
6. Wehler, C. et al.: *The use of solid-phase C1q and horseradish peroxidase-conjugated goat anti-IgG for the detection of immune complexes in human serum.* Molec. Immunol. **18**, 157-162 (1981).
7. Endo, Y. et al.: *An enzyme-linked immunoassay for the measurement of circulating immune complexes.* Ann. Clin. Biochem. **20**, 163-168 (1983).
8. Di Mario U. et al.: *The detection of circulating immune complexes in normal subjects using four different methods.* J Clin Lab Immunol **5** (2), 95-99 (1981).
9. Jans, H. et al.: *Circulating immune complexes in healthy persons.* Scand J Rheumatol **11**(4), 194-198 (1982).

**CIC-C1Q EIA**

**Precoated Microtiter Plate**

Wash 2 x

**100 µl Incubation Buffer,  
Standards, Controls, Plasma  
or Serum Samples (1:50)**

*Incubate 1 hour (± 5 min) at 18-28°C  
on a plate rotator*

Wash 3x

**Add 100 µl Enzyme Label**

*Incubate 30 min (± 5 min) at 18-28°C  
on a plate rotator*

Wash 3x


**Add 100 µl pNPP Substrate**

*Incubate 30 min (± 5 min) at 18-28°C  
On a plate rotator*

**Add 100 µl Stop Solution**

**Read absorbance at 405 nm (within 30 minutes)**

**TIME TO RESULT: 2 HOURS**

**APPENDIX IV**  
**SYMBOLS/ SYMBOLE/ SYMBOLES/ SIMBOLI/ SIMBOLOS**

Symbol	Explanation
	Use By Verwendbar bis Utiliser jusqu'au Utilizzare entro Fecha de caducidad
<b>REF</b>	Catalogue number Bestellnummer Référence du catalogue Numero di catalogo Número de catálogo
<b>LOT</b>	Batch code Chargenbezeichnung Code du lot Codice del lotto Codigo de lote
<b>IVD</b>	<i>In Vitro</i> Diagnostic Medical Device <i>In Vitro</i> Diagnostikum Dispositif médical de diagnostic <i>in vitro</i> Dispositivo medico-diagnostico <i>in vitro</i> Producto sanitario para diagnóstico <i>in vitro</i>
	Contains sufficient for <n> tests Ausreichend für „n“ Ansätze Contenu suffisant pour „n“ tests Contenuto sufficiente per „n“ saggi Contenido suficiente para <n> ensayos
	Consult Instructions for Use- Gebrauchsanweisung beachten Consulter le mode d'emploi Consultare le istruzioni per l'uso Consulte las instrucciones de uso
	Temperature limitation Zulässiger Temperaturbereich Limites de température Limiti di temperatura Limite de temperatura
	Upper limit of temperature Temperaturobergrenze Limite supérieure de température Limite superiore di temperatura Límite superior de temperatura
<b>MP</b>	Microtiterplate Mikrotiter-Platte Microplaque Micropiastra Microplaca

Symbol	Explanation
<b>BUF WASH 10X</b>	Wash Buffer Concentrate (10x) Wasch-Puffer Konzentrat (10x) Concentré de tampon de lavage Tampone di lavaggio concentrato (10x) Tampón de lavado concentrado (x10)
<b>BUF INC</b>	Incubation Buffer Inkubations-Puffer Tampon d'incubation Tampone d'incubazione Tampón de incubación
<b>CALA - CALE</b>	Calibrator A - E Kalibrator A - E Calibreur A - E Calibratore A - E Calibrador A - E
<b>CONTROL L</b>	Control Low Kontrolle tief Contrôle bas Controllo basso Control bajo
<b>CONTROL H</b>	Control High Kontrolle hoch Contrôle élevé Controllo alto Control alto
<b>EL</b>	Enzyme Label Enzym-Marker Marqueur enzymatique Marcato enzimatico Marcador enzimático
<b>SUBS PNPP</b>	pNPP Substrate pNPP-Substrat Substrat pNPP Substrato di pNPP Substrato de pNPP
<b>SOLN STOP</b>	Stop Solution Stopp-Lösung Solution stop Soluzione stoppante Solución de parada

Printing Date  
2012-11-09