# GHB

## Gamma-Hydroxybutyric Acid

### Enzymatic Assay

KK-GHB


Revision date: 2012-05-15

# ENGLISH

## INTENDED USE

The BÜHLMANN gamma-hydroxybutyric acid (GHB) kit has been designed for the direct and quantitative determination of GHB in human urine and serum. The assay can be adapted to clinical chemistry analyzers according to specific instrument protocols.

## PRINCIPLE OF THE ASSAY

The GHB is converted to its metabolite, succinic semialdehyde (SSA), by the action of GHB-Dehydrogenase (Rec. Enzyme) and oxidized nicotinamide adenine dinucleotide (NAD<sup>+</sup>). The increase in absorbance at 340 nm resulting from the reduction of NAD<sup>+</sup> into NADH is proportional to the amount of GHB in the sample.

## REAGENTS SUPPLIED AND PREPARATION

Reagents <sup>1)</sup>	Quantity	Code	Reconstitution
Incubation Buffer	1 vial 12 ml	B-GHB-IB	Ready to use
Cofactor lyophilized NAD <sup>+</sup>	1 vial	B-GHB-CF	Add 5.6 ml of deionized water <sup>4)</sup>
Enzyme lyophilized recombinant Enzyme	2 vials	B-GHB-E	Add 4.2 ml of deionized water <sup>4)</sup> , <b>do not vortex</b>
Calibrators <sup>2)</sup> lyophilized GHB	2 vials	B-GHB-CASET	Add 2 ml of deionized water <sup>4)</sup>
Controls Low / High <sup>3)</sup> lyophilized GHB in human urine	2 vials	B-GHB-CONSET	Add 2 ml of deionized water <sup>4)</sup>

Table 1

- 1) All reagents contain sodium azide (<0.1%) as preservative.
- 2) After reconstitution the GHB concentration of the Calibrators is 10 and 100 mg/L, respectively.
- 3) The Controls contain lot-specific amounts of GHB. Refer to the QC data sheet for actual concentrations.
- 4) For reconstitution see chapter Procedural Notes.

## STORAGE AND SHELF LIFE OF REAGENTS

Unopened Reagents	
The unopened kit components are stable at 2 - 8 °C, <b>except Cofactor that must be stored at -20 °C upon receipt</b> . Do not use reagents after expiration of the kit.	
Opened / Reconstituted Reagents	
Incubation Buffer	At 2-8 °C until the expiration date
Cofactor	Store for up to 2 months at 2-8 °C
Enzyme	
Calibrators	
Controls	

Table 2

## WARNINGS AND PRECAUTIONS

- B-GHB-CASET and B-GHB-CONSET contain GHB at <100 mg/L. At this concentration GHB is considered as non-toxic.
- All reagents contain sodium azide at <0.1%. Sodium azide may react with copper or lead and form explosive compounds in metal drain lines. Upon disposal, flush with a large volume of water to prevent azide build up.
- Do not ingest.
- Use only deionized H<sub>2</sub>O for reconstitution. See chapter Procedural Notes.

- Calibrators and Controls contain components of human urine. Thus the risk of transmitting pathogens of known or unknown origin cannot be completely excluded.
- Urine samples should be handled as potentially infectious.
- All products and samples containing human material should be handled in accordance with good laboratory practice using appropriate precautions.

## MATERIALS REQUIRED BUT NOT PROVIDED

- Accurate pipeting devices (1, 5 or 10 ml pipettes)
- Test tubes/rack
- Clinical chemistry analyzer with an optical filter at 340 nm
- Deionized water
- 0.9 % saline solution

## SPECIMEN COLLECTION AND STORAGE

- Depending on the dead volume of the analyzer, ~110 µl sample volume must be transferred into the instrument specific sample tubes. Less than 10 µl of urine or serum sample is needed to run the test. For details refer to the instrument manual.
- Urine: Collect urine (pH 5-8) without additives, keep the samples refrigerated, and have them filtered or centrifuged. Collect the supernatant.
- Serum: Refer to page 4 to learn about the interference of hemolytic, lipemic or icteric samples. Collect blood into plain tubes (no anti-coagulant), avoid haemolysis, leave to clot for one hour, centrifuge for 10 minutes at approximately 1500 x g at room temperature (18-28 °C), and collect the supernatant.
- Sample storage: Urine and serum samples can be stored at 2-8 °C for at least two weeks. For longer storage keep the samples frozen at -20 °C. GHB stability in samples has intensively been investigated (8).

## PROCEDURAL NOTES

### Reconstitution of Reagents:

- Calibrator and Control vials are under vacuum. Thus, it is mandatory to remove the caps slowly and carefully in order to prevent loss of material.
- Cofactor (B-GHB-CF), Calibrators (B-GHB-CASET) and Controls (B-GHB-CONSET): Vortex the vials for 30 seconds after reconstitution or use a suspension mixer for 15 minutes.
- Enzyme (B-GHB-E): **DO NOT VORTEX!** Use a suspension mixer for 15 minutes at RT until the lyophilized Enzyme has completely dissolved. Alternatively, gently swivel the vial after reconstitution, leave it 15 minutes at RT, and gently move it upside-down afterwards.
- On board, the Enzyme can be kept up to 2 months at a temperature up to 15 °C.

---

## ASSAY PROCEDURE

GHB enzymatic assay is designed for automated analysis on clinical chemistry analyzers. It has to be programmed as a three reagent test. For details refer to the analyzer specific protocols.

The following protocol has been established for KoneLab 30 (Thermo) and is used here as a general example for the test procedure. For applications on other instruments, we will provide assistance.

1. 100 µl of Incubation Buffer (R1)  
+ 8 µl of urine or serum sample (S)  
+ 7 µl of deionized water
2. +50 µl of Cofactor (R2)  
→ Incubate for 1.5-2 min at 37 °C
3. +85 µl of Enzyme (R3)
4. T<sub>0 min.</sub>: Measurement at 340 nm (Blank)  
→ Incubation for 5-6 min at 37 °C
5. T<sub>5-6 min.</sub>: Measurement at 340 nm

---

## RESULTS

### Calibration

Use endpoint mode at 340 nm. The absorbance obtained after reading 1 (M1) and 2 (M2) is calculated using the absorbance differences (M2-M1). The standard curve will be calculated from the two calibrators using a linear regression mode. Refer to Table 11 and Figure 1 for typical results and standard curve. Depending on the instrument it might be necessary to read additionally at a secondary wavelength (e.g. 700 nm).

Calibration is stable for at least 2 weeks. Calibration of the assay should be carried out every 14 days, when a fresh vial of Enzyme is reconstituted or when the controls are out of the confidence range.

### Samples and Controls

The absorbance at 340 nm is recorded for the Calibrators and the GHB concentration for each Control and sample is calculated as indicated above.

---

## CALCULATION

**Automated procedure:** For calculation of the results use the endpoint mode with two calibrators (10 and 100 mg/L). Refer to the instrument manual for further details.

**Dynamic range:** 5 - 230 mg/L

---

## STANDARDIZATION

The GHB Calibrators have been calibrated against commercially available GHB quantified with HPLC.

---

## QUALITY CONTROL

- The values of the Low and High Control provided with the kit must be within the lot specific range indicated on the corresponding QC data sheet. Otherwise, the standard curve has to be re-calibrated.
- It is good laboratory practice to record the following data for each assay: kit lot number, reconstitution dates of kit components, results for Calibrators and Controls, and internal urine or serum pool.
- The precision and expected values of standard curve and controls should be within established limits. The confidence limits for the controls are lot-specific and printed on the QC data sheet delivered with the kit.

- Reliable results will only be obtained by using precise laboratory techniques (current GLP guidelines) and by accurate attention of this instruction for use.

---

## LIMITATIONS

Samples exceeding 230 mg/L should be diluted 1:10 with 0.9% saline solution (1 volume sample + 9 volumes saline solution). The results have to be multiplied by 10.

If the OD of the urine samples exceeds 1.0 OD at 340 nm, results may be falsely positive.

Positive results should be confirmed by alternative methods.

Substances and/or factors other than those investigated in the specificity study may interfere with the test and cause false results.

---

## PERFORMANCE CHARACTERISTICS

Assay performance characteristics have been determined on KoneLab 30.

**Limit of Blank (LoB): <1 mg/L.** The LoB has been established by repeated measurements of blank values (0.9% NaCl, n= 60) in accordance with CLSI protocol EP17-A.

**Limit of Detection (LoD): 1.5 mg/L.** The LoD has been established by repeated measurements of two samples containing 4.8 and 6.2 mg/L GHB (n= 40 and 10, respectively) in accordance with CLSI protocol EP17-A (Table 12).

**Limit of Quantification (LoQ): Urine: 5.0 mg/L; serum: <5.0 mg/L.** The LoQ of urine and serum was determined by three measurements of three urine and serum samples at concentrations between 3 and 16 mg/L. The sensitivity at 10 % CV was determined to be 5.0 and <5.0 mg/L, respectively.

**Precision: Repeatability: ≤5 % CV; Between run: <10 % CV; Between day: <5 % CV; Total precision: <10 % CV.** The precision of urine and serum samples has been determined in accordance with CLSI protocol EP5-A2 by repeated measurements in two runs per day over a period 10 – 20 work days (Table 13).

**Dilution linearity: Urine 100 - 105 %; serum 99 - 100 %.** Three urine and two serum samples with elevated GHB concentration have been diluted with 0.9% NaCl solution (Table 14). The mean value is 103 % in urine and 100 % in serum. It is linear between 5 and 230 mg/l.

**Recovery: Urine: 95 - 107 %; serum: 106 - 113 %.** Three different urine and serum samples have been spiked with increasing amounts of GHB. The mean value is 100 % in urine and 110 % in serum (Table 15).

**Specificity:** The substances listed in Table 16 have been analyzed between 10 and 1000 mg/L to determine the enzyme specificity.

## INTERFERING SUBSTANCES

**Therapeutic drugs and drugs of abuse:** The **therapeutic drugs** and **drugs of abuse** listed in Table 17 and Table 18 were evaluated according to the approved guideline for interference testing in clinical chemistry, EP7-A2 (7), on the KoneLab 30. No interference has been observed up to the listed concentrations.

The interference of **Ethanol** is shown in Figure 2. Up to 3‰ the measured GHB concentration is below 10 mg/L. 1 g/L Ethanol raises the GHB value by 3 mg/L.

**Serum Indices:** No interference is detected with the following substances up to the listed concentrations: **Triglycerides** (Intralipid® 275 mg/dL; equivalent to 7.7 mmol/L triglycerides), **conjugated bilirubin** (360 µmol/L; 30 mg/dL), **unconjugated bilirubin** (513 µmol/L; 30 mg/dL) or **haemoglobin** (3.1 mmol/L; 500 mg/dL) on KoneLab 30.

## EXPECTED VALUES

The following reference range for urine and serum (97.5<sup>th</sup> percentile) has been determined from individuals, who did not consume GHB:

Urine: 97.5 <sup>th</sup> percentile:	3.1 mg/L
Median:	<1.5 mg/L
n =	46
Serum: 97.5 <sup>th</sup> percentile	6.1 mg/L
Median:	<1.5 mg/L
n =	65

## METHOD COMPARISON

34 urine samples have been compared to a published IC-GHB method (6) (Figure 3):

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Controls (2 replicates)	Values (ΔOD)	Replicates [mg/L]	Mean [mg/L]	CV [%]
low	0.017	13.0	13.0	0.3
	0.017	13.1		
high	0.096	76.7	76.7	0.0
	0.095	76.6		

Table 12: LoB , LoD and LoQ

LoB [mg/L]	LoD [mg/L]	LoQ [mg/L]
<1	1.5	Urine: 5 Serum: <5

Table 13: Assay Precision

Sample [mg/L]	Total Precision [%]	Repeatability (Within Run) [%]	Between Run [%]	Between Day [%]
<b>Urine</b>				
11.2	8.6%	4.9%	7.0%	0.0%
55.3	4.5%	3.0%	3.0%	1.5%
100.3	5.2%	1.0%	2.9%	4.2%
<b>Serum</b>				
12.0	9.7%	5.0%	8.3%	0.0%
55.2	4.1%	2.1%	3.6%	0.0%
106.6	3.7%	1.4%	3.5%	0.0%

Table 14: Dilution Linearity

Example 1 of 3	Dilution	Observed [mg/L]	Expected [mg/L]	O/E [%]
Urine 1	1:1	238.6	287.3	83
	1:1.5	170.6	191.5	89
	1:2.25	119.9	127.7	94
	1:3.4	85.1	85.1	100
	1:5.1	58.8	56.7	104
	1:7.6	41.1	37.8	109
	1:11.4	27.8	25.2	110
	1:17.1	19.6	16.8	117
	1:25.6	13.0	11.2	116
	1:38.4	8.3	7.5	112
Mean urine 1				<b>103</b>
Mean urine 2				<b>105</b>
Mean urine 3				<b>100</b>
<b>Mean urine</b>				<b>103</b>
Serum 1	1:1	239.9	295.3	81
	1:1.5	184.0	210.9	87
	1:2.25	132.7	140.6	94
	1:3.4	93.7	93.7	100
	1:5.1	66.1	62.5	106
	1:7.6	45.1	41.7	108
	1:11.4	30.5	27.8	110
	1:17.1	20.0	18.5	108
	1:25.6	12.4	12.3	101
	1:38.4	7.9	8.2	96
1:57.7	5.1	5.5	94	
Mean serum 1				<b>99</b>
Mean serum 2				<b>100</b>
<b>Mean serum</b>				<b>100</b>

Table 15: Spiking Recovery

Example 1 of 3	Native [mg/L]	Spiked with [mg/L]	Observed [mg/L]	Expected [mg/L]	O/E [%]
Urine1	<1.5	10	10.1	10.0	97
		25	15.2	15.0	91
		50	49.1	50.0	98
		100	90.5	100.0	95
Mean urine 1					95
Mean urine 2					98
Mean urine 3					107
<b>Mean urine</b>					<b>100</b>
Serum1	<1.5	10	11.8	10.0	107
		25	18.3	15.0	118
		50	54.7	50.0	110
		100	102.6	100.0	104
Mean serum 1					110
Mean serum 2					113
Mean serum 3					106
<b>Mean serum</b>					<b>110</b>

Table 16: **Enzyme Specificity**

Component/substrate	Substrate specificity
γ-hydroxyvaleric acid, (GHV)	<0.1 %
γ-butyrolactone, (GBL)	4.0 %
1,4-butanediol, (1,4-BD)	<0.1 %
α-hydroxybutyric acid, (AHB)	<0.1 %
β-hydroxybutyric acid, (BHB)	<0.1 %
γ-hydroxyvaleric acid, (GHV)	<0.1 %
succinic acid	<0.1 %
γ-valerolactone, (GVL)	<0.1 %

**Interferences**

Table 17 **Interference of therapeutic drugs**

Therapeutic drugs	No interference up to	
Amikacin	53.2	μmol/L
Caffeine	84.5	μmol/L
Carbamazepine	75.4	μmol/L
Cyclosporine	480	nmol/L
Digoxin	4.43	nmol/L
Ethosuximide	970	μmol/L
Gentamicin	18.6	μmol/L
Lithium	1.77	mmol/L
Lithium (Vitros)	2.25	mmol/L
Methotrexate	7.78	μmol/L
Paracetamol	1.44	mmol/L
Phenobarbitone	223	μmol/L
Phenytoin	74.9	μmol/L
Primidone	58.8	μmol/L
Salicylic acid	2.86	mmol/L
Theophylline	177	μmol/L
Tobramycin	19.5	μmol/L
Valproic acid	998	μmol/L
Vancomycin	20.3	μmol/L

Table 18 **Interference of Drugs of abuse**

Substance	No interference up to	
Amphetamines (d-Methamphetamines)	1250	ng/mL
Barbiturates (Secobarbital)	375	ng/mL
Benzodiazepines (Oxazepam)	390	ng/mL
Cannabinoids (11-Nor-Δ-9-THC-9-COOH)	65	ng/mL
Cocaine (Benzoyllecgonine)	500	ng/mL
LSD	1	ng/mL
Methadone	375	ng/mL
Opiates (Free Morphine)	2500	ng/mL
Phencyclidine	31	ng/mL
Propoxyphene (Norpropoxyphene)	375	ng/mL
Nortriptyline	375	ng/mL

Figure 2

**Interference of Ethanol**



Figure 3 **Correlation between KK-GHB and IC-GHB**

**Table description:** “Calculation” (page 3), “Performance Characteristics” (page 4), and “Interpretation of Results” (page 4).

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1. Bravo DT et al.: *Reliable, sensitive, rapid and quantitative enzyme-based assay for gamma-hydroxybutyric acid (GHB)*. J Forensic Sci. 2004, **49**(2): 379-87
2. Mari F et al.: *What constitutes a normal ante-mortem urine GHB concentration?* J Forensic Leg Med. 2009, **16**(3): 148-51
3. Andresen H et al.: *Liquid ecstasy - a significant drug problem*. Dtsch Arztebl Int. 2008; **105**(36): 599-603
4. Haller C et al.: *GHB urine concentrations after single-dose administration in humans*. J Anal Toxicol. 2006; **30**(6): 360-4
5. Elliott SP: *The many faces of gamma hydroxybutyric acid (GHB)*. Drug Monit. 2003; **4**(2): 9-12
6. Andresen H et al.: *Gamma-hydroxybutyrate in urine and serum: additional data supporting current cut-off recommendations*. Forensic Sci Int. 2010; **200**(1-3): 93-9
7. Hasan L et al.: *An enzymatic method to determine  $\gamma$ -hydroxybutyric acid in serum and urine*. Ther Drug Monit. 2011; **33**(6): 757-65
8. Jordi M et al.: *GHB Determination with Ion Chromatography* Therap. Drug Monit. 25 (2003) 486
9. Clinical Laboratory Standards Institute: *Interference Testing in Clinical Chemistry*; approved guideline (EP7-A2), 2005

Symbol	Explanation	Symbol	Symbol
	Use By Verwendbar bis Utiliser jusqu'au Utilizzare entro Fecha de caducidad	<b>BUFINC</b>	Incubation Buffer Inkubations-Puffer Tampon d'incubation Tampone di incubazione Tampón de incubación
<b>REF</b>	Catalogue number Bestellnummer Référence du catalogue Numero di catalogo Número de catálogo	<b>COF</b>	Cofactor Kofaktor Cofacteur Cofattore Cofactor
<b>LOT</b>	Batch code Chargenbezeichnung Code du lot Codice del lotto Codigo de lote	<b>ENZ</b>	Enzyme Enzym Enzyme Enzime Enzime
<b>IVD</b>	<i>In Vitro</i> Diagnostic Medical Device <i>In Vitro</i> Diagnostikum Dispositif médical de diagnostic <i>in vitro</i> Dispositivo medico-diagnostico <i>in vitro</i> Producto sanitario para diagnóstico <i>in vitro</i>	<b>CAL A</b>	Calibrator A Kalibrator A Calibrateur A Calibratore A Calibrador A
	Contains sufficient for <n> tests Ausreichend für "n" Ansätze Contenu suffisant pour „n“ tests Contenuto sufficiente per „n“ saggi Contenido suficiente para <n> ensayos	<b>CAL B</b>	Calibrator B Kalibrator B Calibrateur B Calibratore B Calibrador B
	Consult Instructions for Use Gebrauchsanweisung beachten Consulter le mode d'emploi Consultare le istruzioni per l'uso Consulte las instrucciones de uso	<b>CONTROL L</b>	Control Low Kontrolle tief Contrôle bas Controllo basso Control bajo
	Temperature limitation Zulässiger Temperaturbereich Limites de température Limiti di temperatura Limite de temperatura	<b>CONTROL H</b>	Control High Kontrolle hoch Contrôle élevé Controllo alto Control alto
	Upper limit of temperature Temperaturobergrenze Limite de température maximale Limite de temperatura maximale Límite de temperatura maximale		

Printing Date  
2012-05-29

Patent pending