Quantum Blue[®] fCAL

Quantitative Lateral Flow Assay

For *In Vitro* Diagnostic Use.

LF-CAL25

25 tests

Release date: 2019-12-20
Version A1

INTENDED USE

The BÜHLMANN Quantum Blue® fCAL is an *in vitro* diagnostic test for the quantitative determination of calprotectin in human stool specimens intended as an aid to the assessment of intestinal mucosal inflammation. The assay results can be used as an aid to diagnosis in distinguishing organic, inflammatory disease of the gastrointestinal tract (inflammatory bowel disease, IBD, specifically Crohn's disease or ulcerative colitis, UC) from functional disease (irritable bowel syndrome, IBS) (ref. 1-7), in patients with chronic abdominal pain.

For laboratory use only.
EU: For professional use.

PRINCIPLE OF THE ASSAY

The test is designed for the selective measurement of calprotectin antigen by sandwich immunoassay. A monoclonal capture antibody (mAb) highly specific for calprotectin is coated onto the test membrane. A second monoclonal detection antibody conjugated to gold colloids is deposited onto the conjugate release pad and released into the reaction system after addition of the extracted and diluted stool sample. The calprotectin/anti-calprotectin gold conjugate binds to the anti-calprotectin antibody coated on the test membrane (test line) and the remaining free anti-calprotectin gold conjugate binds to the goat anti-mouse antibody coated on the test membrane (control line). The signal intensities of the test line (T) and the control line (C) are measured quantitatively by the BÜHLMANN Quantum Blue® Reader.

REAGENTS SUPPLIED AND PREPARATION

Reagents	Quantity	Code	Comments
Test Cassette	25 pieces	B-CAL-TC	vacuum-sealed in a foil pouch
Extraction Buffer	1 bottle 125 mL	B-CAL-EX	Ready to use
Controls Low* / High*	2 vials, 0.5 mL	B-LFCAL-CONSET	Ready to use
RFID Chip Card	1 piece	B-LFCAL-RCC	White plastic card

Table 1

* The controls contain lot specific amounts of native human calprotectin. Refer to the additional QC data sheet for actual concentrations.

STORAGE AND SHELF LIFE OF REAGENTS

All kit components are stable at 2-8 °C until the expiration date printed on the labels.

REAGENTS & MATERIAL SUPPLIED SUPPLEMENTARY

Fecal extraction devices

Fecal extraction devices described below are not delivered with the kit and either of them has to be ordered with the kit.

Extraction device Kits	Quantity	Code
CALEX® Cap device	Packages of 50, 200 or 500 tubes available, filled with 5 mL extraction buffer Ready to use	B-CALEX-C50 B-CALEX-C200 B-CALEX-C500
Smart-Prep	50 tubes consisting of spatulas and base caps	B-CAL-RD
ScheBo® Quick-Prep™	50 tubes consisting of tube, cone & dosing tip, each filled with 1.3 mL extraction buffer Ready to use	B-CAL-SO50

Table 2

MATERIALS REQUIRED BUT NOT PROVIDED

- Vortex mixer for stool extraction
- Precision pipettes with disposable tips: 10-200 µL and 1 mL
- Centrifuge
- 5 mL polypropylene or polystyrene tubes for the dilution of the extracts
- Timer (optional)
- Quantum Blue® Reader available from BÜHLMANN (order code: BI-POCTR-ABS)
- Soft tissues or blotting paper

PRECAUTIONS

Safety precautions

- The controls of this test contain components of human origin. Although tested and found negative for HBV surface antigen, HCV and HIV1/2 antibodies, the reagents should be handled as if capable of transmitting infections and should be handled in accordance with Good Laboratory Practices (GLP) using appropriate precautions.
- Patient specimens should be handled as if capable of transmitting infections and should be handled in accordance with Good Laboratory Practice (GLP) using appropriate precautions.
- Reagents: Avoid contact of reagents with the skin, eyes, or mucous membranes. If contact does occur, immediately wash with generous amounts of water; otherwise, irritation can occur.
- Reagents and chemicals have to be treated as hazardous waste according to the national biohazard safety guideline or regulation.

Technical precautions

Kit components

- All reagents and test samples must be equilibrated to room temperature (18-28 °C) before starting the assay.
- Components must not be used after the expiration date printed on the labels.
- Do not mix different lots of reagents.
- Test cassettes cannot be re-used.

Test procedure

- Read the instructions carefully prior to carrying out the test. Test performance will be adversely affected, if reagents are incorrectly diluted, handled or stored under conditions other than those as detailed in this instruction for use.
- The Quantum Blue® Reader must be switched on and programmed for the Quantum Blue® fCAL assay (CAL_0 or CAL_720) before starting the test (see Quantum Blue® Reader manual).
- Use the white RFID chip card in order to change lot-specific test parameters.
- Patient samples that are not properly handled may cause inaccurate results.
- In order to receive reliable and quantitative results it is important to homogenize the stool sample entirely in the extraction device.
- Diluted samples should be used within several hours and cannot be stored for a longer time period.
- With BÜHLMANN Smart-Prep and ScheBo® Quick-Prep™, it is important to centrifuge the extracts before storage. Centrifuge the tubes for 5 minutes at 3000 x g. After centrifugation the supernatant must be transferred into a fresh storage tube.

SPECIMEN COLLECTION AND STORAGE

For the extraction procedure, less than 1 g of native stool specimen is required. Collect stool specimen into plain tubes.

Important: The specimen must be collected without any chemical or biological additives.

Specimen transport

Stool specimens should be received for processing by the laboratory within 3 days of collection. The specimens may be transported at room temperature or refrigerated.

Specimen storage

Stool specimens should be refrigerated at 2-8 °C and extracted within 3 days of receipt at the laboratory. Do not store samples at elevated temperatures.

Extract storage

Fecal calprotectin extracts obtained with CALEX® Cap device are stable at room temperature for 3 days and at 2-8°C for up to 6 days. For longer storage, freeze extracts at -20°C. Frozen extracts are stable for a period of up to 23 months.

CALEX® Cap extracts can be stored and frozen directly within the CALEX® Cap device. Extracts can be subject to four freeze-thaw cycles. Prior to measurement, allow frozen extracts to equilibrate to room temperature. For re-use / re-

measurement of the extracts see step 2 under the chapter assay procedure.

Fecal calprotectin extracts obtained by manual weighing method, by BÜHLMANN Smart-Prep or by ScheBo® Quick-Prep™ is stable at 2-8 °C for ≤7 days or at -20 °C for 36 months.

ASSAY PROCEDURE

The assay procedure consists of three steps:


1. Extraction of stool samples

The extraction is described in the instruction for use delivered with the respective extraction devices.

CALEX® Cap device: Liquid stool samples can be pipetted directly into the CALEX® Cap device. Unscrew the blue cap and pipet 10 µL of stool sample into the device. Recap the CALEX® Cap device and proceed with vortexing step according to the extraction procedure described and illustrated in the instruction for use delivered with the CALEX® Cap device.

2. Sample processing


- Smart-Prep or ScheBo® Quick Prep™: Let the stool extract settle for 10 minutes after extraction. Dilute the supernatant 1:16 with extraction buffer (e.g. 20 µL extract and 300 µL extraction buffer) and mix well. Let the samples equilibrate for at least 5 minutes to 18-28 °C prior to proceeding to the next step (step no. 3).
- CALEX® Cap device: After extraction, let the stool extract settle for 10 minutes with the white head of the device down. Unscrew the blue cap and dilute the supernatant 1:1.6 with extraction buffer (e.g. 200 µL extract and 120 µL extraction buffer) and mix well. Let the samples equilibrate for at least 5 minutes to 18-28 °C prior to proceeding to the next step (step no. 3).

3. Lateral flow assay procedure and readout


There are two alternative methods available on the Quantum Blue® Reader: <CAL_720> and <CAL_0>. Select one of these methods before starting the experiments.

Load the lot specific parameters from the white RFID chip card.


3.1. Method <CAL_720> with internal timer

- Load the test cassette onto the test cassette holder of the Reader.
- Add 60 µL of diluted stool extract onto the sample loading port of the test cassette.
- Close the cassette holder and start the measurement by pressing the start button.
- The scan starts automatically after 12 minutes (720 seconds).
- For low / high controls: Repeat step 3.1 using 60 µL of controls instead of diluted stool extract.

3.2. Method <CAL_0> without internal timer

- Add 60 µL of diluted stool extract onto the sample loading port of the test cassette.
- Incubate for 12 minutes +/- 1 minute (set a timer manually).
- Load the test cassette onto the test cassette holder of the Reader.
- Scan the cassette with the Quantum Blue® Reader by pressing the start (<ENTER>) button immediately.
- For low / high controls: Repeat step 3.2 using 60 µL of control instead of diluted stool extract.

Remark: Please refer to the Quantum Blue® Reader manual to learn more about the basic functions and how to initialize and operate the reader, especially how to select test methods, and how to load lot specific parameters from the RFID chip card on the Quantum Blue® Reader.

QUALITY CONTROL

- If the precision of the assay does not correlate with the established limits and repetition excludes errors in technique, check the following issues: i) pipetting,

temperature controlling and timing ii) expiration dates of reagents and iii) storage and incubation conditions.

- Result of the self-test of the Quantum Blue® Reader performed at startup has to be valid.

VALIDATION OF RESULTS

- For a valid test result, the control line (C) must be visible in any case (see figures 1A and 1B). It is used as functional test control only and cannot be used for the interpretation of the test line (T). If the test line (T) is not detectable after 12 minutes of incubation time (figure 1A), the concentration of calprotectin present in the stool sample is below the detection limit. If a test line (T) is detectable after 12 minutes of incubation time (figure 1B), the calprotectin concentration present in the stool sample is calculated by the Quantum Blue® Reader.
- If only the test line (T) is detectable after 12 minutes of incubation time (figure 1C), the test result is invalid and the calprotectin assay has to be repeated using another test cassette.
- If neither the control line (C) nor the test line (T) are detectable after 12 minutes of incubation time (figure 1D), the test result is invalid and the calprotectin assay has to be repeated using another test cassette.
- As the Quantum Blue® Reader allows a quantitative evaluation of the test (T) and control (C) lines, an additional validity check of the control line (C) is undertaken. If the signal intensity of the control line (C) is below a threshold after 12 minutes of incubation time, the test result is also invalid and the Quantum Blue® fCAL assay has to be repeated using another test cassette.

STANDARDIZATION

- The Quantum Blue® fCAL is standardized against the BÜHLMANN fCAL® ELISA (order code: EK-CAL).
- The Quantum Blue® Reader uses a lot-specific standard curve to calculate the calprotectin concentration. This lot-specific standard curve is generated with the median values ($n \geq 10$ measurements each) from at least 10 calibration points obtained from different stool samples with known calprotectin concentrations. The assay range is between 30 and 300 µg/g.
- For quantitative measurements, unknown samples reading above 300 µg/g can be re-tested in the BÜHLMANN Quantum Blue® fCAL high range (order code: LF-CHR25).

LIMITATIONS

- Reagents delivered with the BÜHLMANN Quantum Blue® fCAL kit are intended for the determination of calprotectin levels in human stool samples only.
- Fecal calprotectin values are intended as an aid to diagnosis in distinguishing organic disease from functional disease. Results should always be interpreted in combination with other clinical and laboratory findings.
- Patients who are taking NSAIDs regularly may have elevations in their fecal calprotectin levels.

- Results may not be clinically applicable to children less than 4 years of age who have mildly increased fecal calprotectin levels (ref. 8-11).

INTERPRETATION OF RESULTS

DISTINGUISHING ORGANIC DISEASE FROM FUNCTIONAL GASTROINTESTINAL DISEASE

The determination of fecal calprotectin levels can be used as a reliable and simple aid in distinguishing organic from functional gastrointestinal diseases (ref. 1-7).

The result categories are based on data from clinical studies performed by BÜHLMANN and are BÜHLMANN's recommendations. All test results should be interpreted in conjunction with information available from the patient's clinical symptoms, medical history, and other clinical and laboratory findings.

Clinical thresholds

The following data were established with the BÜHLMANN fCAL® ELISA (order code: EK-CAL).

Results from 58 clinical samples from patients diagnosed with IBS and 131 clinical samples from patients diagnosed with IBD, from an international clinical study, were analyzed to obtain the values described in table 3.

Calprotectin concentration	Interpretation	Follow-up
< 80 µg/g	Normal	None
80 - 160 µg/g	Gray-zone/Borderline	Follow-up within 4 – 6 weeks
> 160 µg/g	Elevated	Repeat as needed

Table 3

Calprotectin values below 80 µg/g

Fecal calprotectin values <80 µg/g are not indicative of inflammation in the gastrointestinal tract. Patients with low calprotectin levels are not likely to be in need of invasive procedures to determine the inflammation cause.

Calprotectin values between and equal to 80 and 160 µg/g

Mid-fecal calprotectin levels between and equal to 80 and 160 µg/g, also called gray-zone levels, are not directly indicative of an active inflammation requiring immediate follow-up with invasive testing. However, the presence of inflammation cannot be excluded. Re-evaluation of fecal calprotectin levels after 4 to 6 weeks is recommended to determine the inflammatory status.

Calprotectin values above 160 µg/g

Fecal calprotectin values >160 µg/g are indicative of neutrophil infiltrate in the gastrointestinal tract; therefore, this may signal the presence of active inflammatory disease. Appropriate further investigative procedures by specialists are suggested to achieve an overall clinical diagnosis.

Clinical evaluation

The ability of the BÜHLMANN fCAL® ELISA to discriminate between patients with IBD and other non-inflammatory GI disorders, including IBS, was tested in a clinical study with a total of 337 adult and pediatric patients. One hundred and thirty five (135) patients had a final diagnosis of IBD (Crohn's disease, ulcerative colitis or indeterminate colitis), 130 patients suffered from IBS and 72 patients presented

with abdominal pain and/or diarrhea, or other GI-related non-inflammatory conditions (refer to table 4). Final diagnosis was supported by endoscopic as well as other clinical findings.

A clinical sensitivity of 93.3% at 80 µg/g and a clinical specificity of 83.7% at 160 µg/g, can be reached in the differentiation between IBD and GI-related non-inflammatory conditions, including IBS. ROC curve analysis resulted in an AUC of 0.923 (refer to table 5).

A clinical sensitivity of 93.3% at 80 µg/g and a clinical specificity of 85.4% at 160 µg/g, can be reached in the differentiation between IBD and IBS. ROC curve analysis resulted in an AUC of 0.933 (refer to table 6).

The optimal cut-off combination for these patient pools could be defined by ROC analysis at 80 µg/g and 160 µg/g calprotectin, which is slightly more stringent than a combination of a **more sensitive lower cut-off of 50 µg/g** with lower performance in specificity, and an **upper cut-off of 200 µg/g** with slightly lower sensitivity (table 7 and 8).

PERFORMANCE CHARACTERISTICS

Limit of Blank (LoB): <7 µg/g calprotectin

The LoB has been established in accordance with CLSI protocol EP17-A in three independent runs using three different lots of test cassettes with 60 blank values in total by using extraction buffer as a sample.

Limit of Detection (LoD): 15 µg/g calprotectin

The LoD has been established with two stool extracts with concentrations of 14.5 and 28.5 µg/g. The samples were measured with three different lots of test cassettes in three independent runs of 10 replicates each, set up within 10 minutes. The averaged SD values were determined and the LoD has there from been calculated in accordance with CLSI protocol EP17-A.

Limit of Quantification (LoQ):

Lower LoQ: ≤30 µg/g calprotectin

Upper LoQ: ≥300 µg/g calprotectin

The LoQ has been established with seven stool extracts at concentrations between 14.5 and 622 µg/g. The samples were measured with three different lots of test cassettes in three independent runs of 10-20 replicates each, set up within 10 minutes and averaged. The resulting precision profile is shown in figure 2. The limit of quantification corresponds to the concentration of calprotectin with an imprecision below 25% CV allowing a quantitative measurement within the range from 30 (lower LoQ) to 300 µg/g (upper LoQ).

Linearity

Four stool samples with elevated calprotectin concentrations were extracted according to the assay procedure. The extracts were diluted with extraction buffer and each dilution was subsequently assayed on two test cassette according to the assay procedure and averaged for each test point. The results showed linearity within the indicated measuring range of 30 to 300 µg/g of the Quantum Blue® fCAL assay for all 4 samples ($R^2 = 0.986-0.995$). One example is shown in (figure 3). The average difference between measured and expected concentrations was 2% corresponding to a recovery of 98%. The recovery varied between 84-114%.

Recovery: 92%

Four stool extracts containing low calprotectin concentrations were spiked with different amounts from a highly concentrated calprotectin stock solution. Each sample was measured before and after spiking on three test cassettes according to the assay procedure. The three replicates were averaged and the mean values are presented in table 9. The recovery varied 83-105%.

Precision (repeatability): 15.5-22.2% CV

The repeatability of the Quantum Blue® fCAL assay was calculated from 2 extracted stool samples containing 49 and 138 µg/g calprotectin. Each sample was tested according to the assay procedure in one run of 20 replicates set up within 10 minutes. The mean values of three different lots of test cassette are presented in table 10. The repeatability varied between 15.5-22.2% CV.

Inter-lot Precision: 20.2-24.5% CV

The inter-lot precision of the Quantum Blue® fCAL assay was calculated from 4 extracted stool samples containing between 30 and 244 µg/g calprotectin. Each sample was tested according to the assay procedure in three independent runs of 10-20 replicates each, set up within 10 minutes using three different lots of test cassette. The inter-lot data of the 4 samples varied from 20.2 to 24.5% CV.

Method Comparison: $R^2 = 0.806$; $y = 1.02x - 11.1$ µg/g

One hundred and one (101) stool extracts from subjects undergoing IBD vs. IBS screening exhibiting calprotectin concentrations within the indicated measuring range of the Quantum Blue® fCAL assay were analyzed according to the assay procedure using a data set originated from test cassette of six different lots and compared with the values obtained by the BÜHLMANN fCAL® ELISA (order code: EK-CAL). The correlation data are illustrated in figure 4.

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Figure 1D

Figure 1

Clinical study – distinguishing organic disease from functional gastrointestinal disease

Final diagnosis	Distribution of patients' results in numbers (percent) within BÜHLMANN fCAL® ELISA diagnostic ranges.			
	< 80 µg/g	80 - 160 µg/g	> 160 µg/g	Total
IBD	9 (6.7%)	12 (8.9%)	114 (84.4%)	135 (100%)
IBS	94 (72.3%)	17 (13.1%)	19 (14.6%)	130 (100%)
Other GI	48 (66.7%)	10 (13.9%)	14 (19.4%)	72 (100%)

Table 4

IBD vs. non-IBD	Clinical decision point	
	80 µg/g	160 µg/g
Sensitivity (95% CI)	93.3% (87.7%, 96.9%)	84.4% (77.2%, 90.1%)
Specificity (95% CI)	70.3% (63.5%, 76.5%)	83.7% (77.8%, 88.5%)
PPV (95% CI)	67.7% (60.5%, 74.4%)	77.6% (69.9%, 84.0%)
NPV (95% CI)	94.0% (89.0%, 97.2%)	88.9% (83.6%, 93.0%)
ROC AUC (95% CI)	0.923 (0.893, 0.953)	

Table 5

IBD vs. IBS	Clinical decision point	
	80 µg/g	160 µg/g
Sensitivity (95% CI)	93.3% (87.7%, 96.9%)	84.4% (77.2%, 90.1%)
Specificity (95% CI)	72.3% (63.8%, 79.8%)	85.4% (78.1%, 91.0%)
PPV (95% CI)	77.8% (70.6%, 83.9%)	85.7% (78.6%, 91.2%)
NPV (95% CI)	91.3% (84.1%, 95.9%)	84.1% (76.7%, 89.9%)
ROC AUC (95% CI)	0.933 (0.902, 0.963)	

Table 6

IBD vs. non-IBD	Clinical decision point	
	50 µg/g	200 µg/g
Sensitivity (95% CI)	96.3% (91.6%, 98.8%)	80.7% (73.1%, 87.0%)
Specificity (95% CI)	59.9% (52.8%, 66.7%)	87.1% (81.7%, 91.4%)
PPV (95% CI)	61.6% (54.7%, 68.2%)	80.7% (73.1%, 87.0%)
NPV (95% CI)	96.0% (91.0%, 98.7%)	87.1% (81.7%, 91.4%)

Table 7

IBD vs. IBS	Clinical decision point	
	50 µg/g	200 µg/g
Sensitivity (95% CI)	96.3% (91.6%, 98.8%)	80.7% (73.1%, 87.0%)
Specificity (95% CI)	59.2% (50.3%, 67.8%)	90.0% (83.5%, 94.6%)
PPV (95% CI)	71.0% (63.9%, 77.5%)	89.3% (82.5%, 94.2%)
NPV (95% CI)	93.9% (86.3%, 98.0%)	81.8% (74.5%, 87.8%)

Table 8

Non-IBD – IBS + other GI

CI – confidence interval

PPV – positive predictive value

NPV – negative predictive value

ROC AUC – area under receiver operating characteristic curve

Precision profile


Figure 2

Linearity plot


Figure 3

APPENDIX I

TABLES AND FIGURES / TABELLEN UND ABBILDUNGEN / TABLES ET FIGURES / TABELLE E FIGURE / TABLAS Y FIGURAS/ TABELAS E FIGURAS

Spiking recovery

Sample	µg/g	spiked with [µg/g]	Observed [µg/g]	Expected [µg/g]	O/E [%]
P1	53	25	62	78	79
		50	154	103	150
		100	143	153	94
		200	239	253	95
	Mean				105
P2	81	25	107	106	101
		50	120	131	91
		100	162	181	89
		200	265	281	94
	Mean				94
P3	56	25	77	81	94
		50	120	106	113
		100	110	156	71
		200	186	256	73
	Mean				88
P4	160	25	153	185	82
		50	191	210	91
		100	190	260	73
		200	305	360	85
	Mean				83
Mean					92

Table 9

Precision (repeatability)

Lot	Sample	Mean [µg/g]	SD [µg/g]	CV [%]
A	S1	47.8	10.6	22.2
A	S2	141.3	22.0	15.5
B	S1	50.1	10.0	20.1
B	S2	145.5	23.5	16.2
C	S1	50.0	10.3	20.6
C	S2	146.6	31.5	21.7
Mean				19.4

Table 10

Method comparison


Figure 4

APPENDIX II

REFERENCES / LITERATURREFERENZEN / REFERENCES / RIFERIMENTI / REFERENCIAS

1. Fagerhol MK: *Calprotectin, a faecal marker of organic gastrointestinal abnormality*. Lancet 356, 1783-4 (2000)
2. Tibble JA et al.: *A simple method for assessing intestinal inflammation in Crohn's disease*. Gut 47, 506-513 (2000).
3. Tibble JA et al.: *Use of surrogate markers of inflammation and Rome criteria to distinguish organic from nonorganic intestinal disease*. Gastroenterol 123, 450-460 (2002)
4. Jahnsen J, Røseth AG, Aadland E.: *Measurement of calprotectin in faeces..* Tidsskr Nor Lægeforen 128, 743-5 (2008)
5. Manz M et al.: *Value of fecal calprotectin in the evaluation of patients with abdominal discomfort: an observational study*. BMC Gastroenterology 12, 5 (2012).
6. Pavlidis P. et al.: *Diagnostic accuracy and clinical application of faecal calprotectin in adult patients presenting with gastrointestinal symptoms in primary care*. Scand J Gastroenterol. 48, 1048-54 (2013)
7. Konikoff MR and Denson LA: *Role of fecal calprotectin as a biomarker of intestinal inflammation in inflammatory bowel disease*. Inflamm Bowel Dis 12(6), 524-34 (2006)
8. Fagerberg UL et al.: *Colorectal inflammation is well predicted by fecal calprotectin in children with gastrointestinal symptoms*. J Pediatr Gastroenterol Nutr 40, 450-5 (2005)
9. Li F. et al.: *Fecal Calprotectin Concentrations in Healthy Children Aged 1-18 Months*. PLoS ONE 10(3) (2015)
10. Zhu Q. et al. *Fecal Calprotectin in Healthy Children Aged 1-4 Years*. PLoS ONE 11 (3) (2016)
11. Peura S. et al.: *Normal values for calprotectin in stool samples of infants from the population-based longitudinal born into life study*. Scandinavian Journal of Clinical and Laboratory Investigation 78(1-2), 120-124 (2018)

APPENDIX III

INCIDENT REPORTING IN EU MEMBER STATES

If any serious incident in relation to this device has occurred, please report without delay to the manufacturer and competent authority of your Member State.

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