

# BÜHLMANN GanglioCombi® Light ELISA

## Pre-Analytics

<b>Specimen</b>	Serum
<b>Specimen collection</b>	Use conventional plain tubes to collect blood. Process to serum according to manufacturer's instruction.
<b>Storage of Serum</b>	<ul style="list-style-type: none"> <li>• Prepare aliquots</li> <li>• 2-8°C (up to 16 days)</li> <li>• -20°C (up to 12 months)</li> <li>• Freeze/thaw cycles not recommended</li> </ul>

## Number of Analyses


- Kit allows for full flexibility regarding:
  - number of samples to be analyzed
  - use and arrangement of enzyme conjugates on microtiter plates
- Generally, a sample can be investigated for 3 different antibodies with IgG and IgM and/or with IgG/IgM mix conjugate.
- Each kit contains 1 microtiter plate.

Enzyme Conjugate	Number of Analyses per Kit	Number of Analyses per individual Sample
IgG & IgM	12 + 12	2 x 3 (2 conjugates/3 gangliosides)
IgG/IgM Mix	24	1 x 3 (1 mix conjugate/3 gangliosides)

## Literature References

Kuijf M et al., Journal of the Neurological Sciences, 2005  
 Spatola M et al., Neurology, 2016  
 Franciotta D et al., Clin Chem Lab Med, 2018

## Microtiter Plate Layout



- Rows coated with 3 different gangliosides
- Break-away strip convenience

## Reagent Preparation

Reagent	Temperature °C
Wash buffer	2-8
Incubation Buffer	2-8
Enzyme conjugate	2-8
TMB substrate	18-28
Stop solution	18-28

## Sample Preparation

1. Dilute serum samples 1:50 with incubation buffer and vortex gently.
2. Leave diluted samples, reconstituted calibrators and controls at 2-8°C for 30 minutes prior pipetting (for targeted investigation).

## ELISA Procedure\*

<p><b>Precoated Microtiter Plate</b></p> <p>↓ ↻ wash 2x (≥ 300 µL wash buffer)</p> <p><b>Calibrators, Controls, Serum (1:50) (100 µL)</b></p> <p>↓ ↻ incubate (2 hr ± 5 min)</p> <p>↓ ↻ wash 3x (≥ 300 µL wash buffer)</p> <p><b>Enzyme Label (100 µL)</b></p> <p>↓ ↻ incubate (2 hr ± 5 min)</p> <p>↓ ↻ wash 3x (≥ 300 µL wash buffer)</p> <p><b>TMB Substrate (100 µL)</b></p> <p>↓ ↻ incubate (30 min ± 2 min) plate shaker: (400-600 rpm)</p> <p><b>Stop Solution (100 µL)</b> ≤ 30 min</p> <p>↓</p> <p><b>Absorbance (450 nm)</b></p>	(2 - 8°C)
	(18 - 28°C)

\*Please see the full instructions for use (IFU) on [buhlmannlabs.com](http://buhlmannlabs.com) for more detailed information.


BÜHLMANN GanglioCombi® Light ELISA. Not for use in diagnostic procedures. For product availability outside the US, please contact BÜHLMANN Laboratories for further information.

# BÜHLMANN GanglioCombi® Light ELISA

## Methodology


<b>Method</b>	ELISA for Research Use Only
<b>Analyte</b>	anti-ganglioside antibodies: GM1, GD1b, GQ1b
<b>Test</b>	multi-parametric, semi-quantitative
<b>Suggested Use</b>	Targeted investigation of anti-ganglioside antibodies with IgG, IgM, and/or IgM/IgG Mix.

## Principle of the Assay



\*We use double-quality control-checked gangliosides for the coating of the microtiter plates.

## Reference Interval



Healthy Individuals: N = 120 (m = 60; f = 60)  
Negative: 95.0%; Grey zone: 3.9%; Positive: 1.1%

## Performance Data

<b>Crossreactivity:</b>	No (incl. AIDS <sup>1</sup> and ONDs <sup>2</sup> )
<b>Reproducibility (%CV):</b>	7.7-19.1%
<b>Specificity (95 %CI):</b>	79.8% (66.8-93.2%)
<b>Sensitivity (95 %CI):</b>	57.0% (38.6-75.4%)

<sup>1</sup> Autoimmune diseases (N = 50): Vasculitis and ANCA positive denoted: 3 and 10; SLE: 5; RA: 9 Sjögren Syndrome: 6; Others (ANA positive denoted): 3; Autoimmune Thyroiditis: 5; MCTD: 1; Autoimmune peripheral neuropathies (GM1 and GQ1b pos): 1; MG: 7

<sup>2</sup> Other neurological diseases (N = 34): Alcoholic: 1; diabetic: 5; ALS: 15; Sarcoidosis: 4; WM: 4; Chagas Disease: 5

## Results

- The measurement of absorbance is proportional to the titer of ganglioside antibodies in a given sample.
- Titers are expressed as a % ratios of the calibrator

## Handling of Titer levels

Titer levels [%Ratio]	Interpretation [titer]
< 30	negative
30 - 50	grey zone (equivocal)
> 50	positive

## Ordering information

•BÜHLMANN GanglioCombi® Light ELISA 3 gangliosides EK-GCL-S-U

### Variants:

•BÜHLMANN GanglioCombi® MAG ELISA 5 gangliosides, 1xMAG EK-GCM-U

•anti-GM1 Antibodies ELISA 1 ganglioside EK-GM1-GM -U

BÜHLMANN GanglioCombi® Light ELISA. Not for use in diagnostic procedures. For product availability outside the US, please contact BÜHLMANN Laboratories for further information.

CE 0123

**BÜHLMANN**

BÜHLMANN Laboratories AG  
Switzerland  
info@buhlmannlabs.ch  
www.buhlmannlabs.ch  
PHone +41 61 487 12 12

BÜHLMANN Diagnostics Corp  
USA  
info@buhlmannlabs.com  
www.buhlmannlabs.com  
Phone 844 300 9799