

Quantum Blue® Adalimumab

trough level rapid test

Adalimumab
rapid test



- First quantitative rapid test for Adalimumab measurement
- Time to result 15 minutes
- Measuring Range 1.3-35 µg/mL
- High agreement to commercially available ELISAs

For Research Use Only in the US

Quantum Blue® Adalimumab is available for Research Use Only in the US. Not for use in diagnostic procedures.


Quantum Blue® Adalimumab rapid test

Key features

- The first quantitative rapid test for adalimumab
- Measuring range 1.3 -35 µg/mL
- Covers the range of the relevant 5 to 12 µg/mL for trough levels*
- High correlation and agreement to commercially available ELISA assays

* Mitrev et al. 2017

Correlations to ELISA assay:



— Passing-Bablok fit

		Commercially available ELISA test (µg/mL)		
		low	optimal	high
Quantum Blue® Adalimumab (µg/mL)	low	20.0%	5.0%	0.0%
	optimal	0.0%	27.5%	0.0%
	high	0.0%	12.5%	35.0%

total agreement 82.5%



BÜHLMANN Laboratories AG
Switzerland
info@buhlmannlabs.ch
www.buhlmannlabs.ch
Phone +41 61 487 12 12

BÜHLMANN Diagnostics Corp
USA / Canada
info@buhlmannlabs.com
www.buhlmannlabs.com
Phone 844 300 9799

Quantum Blue® Adalimumab (LF-TLAD25-U/ LF-TLAD10-U) is available For Research Use Only in the US. Not for use in diagnostic procedures.

Quantum Blue® is a registered trademark of BÜHLMANN in many countries.

Ordering codes:
LF-TLAD10-U 10 cassettes
LF-TLAD25-U 25 cassettes

Related products:
LF-TLIF10-U 10 cassettes
LF-TLIF25-U 25 cassettes
LF-CAL25-U 25 cassettes
LF-CALE25-U 25 cassettes
LF-CHR25-U 25 cassettes