

Quantum Blue[®] Infliximab trough level rapid test

IFX rapid test
act immediately
before next
infusion



- First rapid test for IFX monitoring
- Time to result 15 minutes
- Measuring range
0.4 to 20 µg/ml IFX in serum
- Highly specific
monoclonal antibody


Health Canada Licence:98838

Quantum Blue® Infliximab rapid test

Key features

- The only quantitative rapid test for infliximab measurements
- Measuring range 0.4 - 20 µg/mL
- Covers the decisive range of 3 to 7 µg/mL for trough levels
- Incorporated in clinical trials in IBD and rheumatoid arthritis
- Highly specific monoclonal antibody
- Very good correlation to established third party ELISA assays

Trough levels infliximab



*Vande Castele et al. 2015



BÜHLMANN Laboratories AG
Switzerland
info@buhlmannlabs.ch
www.buhlmannlabs.ch
Phone +41 61 487 12 12

BÜHLMANN Diagnostics Corp
USA / Canada
info@buhlmannlabs.com
www.buhlmannlabs.com
Phone 844 300 9799

*Health Canada Licence: 98838
Quantum Blue® Infliximab (LF-TLIF25)

Ordering code:
Quantum Blue® Infliximab LF-TLIF25 25 cassettes

Related products:
LF-CAL25 25 cassettes
LF-CALE25-U 25 cassettes
LF-CHR25-U 25 cassettes