

Quantum Blue® Infliximab

trough level rapid test

IFX rapid test



- First rapid test for IFX measurement
- Time to result 15 minutes
- Measuring range
0.4 to 20 µg/mL IFX in serum
- Extended range for measurements during induction
- Compatible with Biosimilars

For Research Use Only in the US

Quantum Blue® Infliximab is available for Research Use Only in the US. Not for use in diagnostic procedures.

Quantum Blue® Infliximab rapid test

Key features

- Quantitative rapid test for infliximab measurements
- Measuring range 0.4 -20 µg/mL, extendable for the study of induction phase
- Compatible with Remsima®, Inflectra® and Flixabi®
- Very good correlation to established third party ELISA assays



BÜHLMANN Laboratories AG
Switzerland
info@buhlmannlabs.ch
www.buhlmannlabs.ch
Phone +41 61 487 12 12

BUHLMANN Diagnostics Corp
USA / Canada
info@buhlmannlabs.com
www.buhlmannlabs.com
Phone 844 300 9799

Quantum Blue® Infliximab (LF-TLIF10-U/
LF-TLIF25-U) is available For Research Use
Only in the US. Not for use in diagnostic
procedures.

Quantum Blue® is a registered trademark of
BÜHLMANN in many countries.

Ordering code:

LF-TLIF10-U	10 cassettes
LF-TLIF25-U	25 cassettes

Related products:

LF-TLAD10-U	10 cassettes
LF-TLAD25-U	25 cassettes
LF-CAL25-U	25 cassettes
LF-CALE25-U	25 cassettes
LF-CHR25-U	25 cassettes