

Quantum Blue® Infliximab

trough level rapid test

IFX rapid test
act immediately
before next
infusion



- First rapid test for IFX monitoring
- Time to result 15 minutes
- Measuring range
0.4 to 20 µg/mL IFX in serum
- Extended range for measurements during induction
- Validated for biosimilars


Health Canada Licence: 98838

Quantum Blue® Infliximab rapid test

Key features

- Quantitative rapid test for infliximab measurements
- Measuring range 0.4 -20 µg/mL, extendable for induction phase
- Covers the range of 3 to 7 µg/mL for trough levels
- Incorporated in clinical trials in IBD and rheumatoid arthritis
- Validated for Remsima®, Inflectra® and Flixabi®
- Very good correlation to established third party ELISA assays

Trough levels infliximab



* Vande Casteele et al. 2015



BÜHLMANN Laboratories AG
Switzerland
info@buhlmannlabs.ch
www.buhlmannlabs.ch
Phone +41 61 487 12 12

BÜHLMANN Diagnostics Corp
USA / Canada
info@buhlmannlabs.com
www.buhlmannlabs.com
Phone 844 300 9799

Health Canada Licence: 98838
Quantum Blue® Infliximab (LF-TLIF10/LF-TLIF25)

Quantum Blue® is a registered trademark of BÜHLMANN in many countries.

Ordering codes:
LF-TLIF10 10 cassettes
LF-TLIF25 25 cassettes
Related products:
LF-TLAD10 10 cassettes
LF-TLAD25 25 cassettes
LF-CAL25 25 cassettes
LF-CAL25 25 cassettes
LF-CHR25 25 cassettes