

# CALEX<sup>®</sup> Cap

## Fecal Extraction Device

Easy to use  
Calprotectin  
Extraction



**NEW**

- Redesigned for total laboratory automation (TLA) systems
- Improved ergonomometry for manual handling
- Fast and easy stool extraction
- High correlation to gold standard manual extraction
- Optimal dilution offers excellent efficiency in extraction

**NEW**

**US:** Available for Research Use Only.  
Not for use in diagnostic procedures.


# CALEX® Cap Unique Features

## New Design




Redesigned sampling pinhead enables total laboratory automation and smooth rack filling

Improved cap ergonometry to facilitate automated and manual decapping



Comparison between the redesigned Calex® Cap and the previous version.  $r^2 = 0.96$



Comparison of Calex® Cap and the gold standard manual extraction.  $r^2 = 0.9$

### SIMPLICITY

Improved ergonometry further supports the simple and safe use of Calex® Cap, significantly reducing hands-on time

### AUTOMATION

CALEX® Cap is redesigned for the use on total laboratory automation (TLA) solutions with full compatibility for the BÜHLMANN fCAL® turbo assay

### RELIABILITY

Excellent precision in stool sampling and high correlation to the gold standard manual stool extraction procedure

### STABILITY

Extensive validation proves stability at RT for 3 days. The Calex® Cap is suitable for air transportation according to IATA 650 (UN3373)

### EFFICIENCY

Optimal extract dilution of 1:500

### EXCLUSIVITY

CALEX® Cap is designed for all BÜHLMANN calprotectin assays



BÜHLMANN Laboratories AG  
Switzerland  
info@buhlmannlabs.ch  
www.buhlmannlabs.ch  
Phone +41 61 487 12 12

BÜHLMANN Diagnostics Corp  
USA / Canada  
info@buhlmannlabs.com  
www.buhlmannlabs.com  
Phone 844 300 9799

CALEX Cap® (B-CALEX-C50-U, B-CALEX-C200-U, B-CALEX-C500-U) is available For Research Use Only in the US. Not for use in diagnostic procedures.

#### Ordering Codes:

B-CALEX-C50-U	50 devices
B-CALEX-C200-U	200 devices
B-CALEX-C500-U	500 devices

#### Related Product:

B-CAL-SAMPLER 25	1 kit
------------------	-------

BÜHLMANN fCAL® and Calex® are registered trademarks of BÜHLMANN in many countries.