

# Quantum Blue® Adalimumab trough level rapid test

## Adalimumab Rapid Test

- First rapid test for Adalimumab measurement
- Time to result 15 minutes
- Measuring Range 1.3-35 µg/mL
- High agreement to commercially available ELISAs

### For Research Use Only in the US and Canada

Quantum Blue® Adalimumab is available for Research Use Only. Not for use in diagnostic procedures.


# Quantum Blue® Adalimumab rapid test

## Key features

- The only quantitative rapid test for adalimumab
- Measuring range 1.3 - 35 µg/mL
- Covers the range of the relevant 5 to 12 µg/mL for trough levels\*
- High correlation and agreement to commercially available ELISA assays

\* Mitrev et al. 2017

## Correlations to ELISA assay:



		Commercially available ELISA test (µg/mL)		
		low	optimal	high
Quantum Blue® Adalimumab (µg/mL)	low	20.0%	5.0%	0.0%
	optimal	0.0%	27.5%	0.0%
	high	0.0%	12.5%	35.0%

— Passing-Bablok fit

total agreement 82.5%



BÜHLMANN Laboratories AG  
Switzerland  
info@buhlmannlabs.ch  
www.buhlmannlabs.ch  
Phone +41 61 487 12 12

BÜHLMANN Diagnostics Corp  
USA / Canada  
info@buhlmannlabs.com  
www.buhlmannlabs.com  
Phone 844 300 9799

Quantum Blue® Adalimumab (LF-TLAD25-U) is available For Research Use Only in the US and Canada. Not for use in diagnostic procedures.

**Ordering code:**  
LF-TLAD25-U 25 cassettes

**Related products:**  
LF-TLIF25-U 25 cassettes  
LF-CALE25-U 25 cassettes